



# cubeX<sup>2</sup>

## cubic zirconia

**TRANSLUCENT LIKE LITHIUM DISILICATE • STRENGTH OF ZIRCONIA**




**TRY IT FOR FREE**  
Ask us for details

- ✓ Beauty and strength (720 MPa)
- ✓ Conservative preparation - 0.8 mm chamfer or rounded shoulder
- ✓ Conventional cementation
- ✓ Warm color match (no graying)
- ✓ Single crowns and 3-unit bridgework in all tooth positions

Available for purchase in 98.5mm discs:

Size	Price	Approximate unit cost (based on 20 units per disc)
98.5mm x 10mm	\$200	\$10.00
98.5mm x 12 mm	\$217	\$10.85
98.5mm x 14mm	\$233	\$11.65
98.5mm x 16mm	\$250	\$12.50
98.5mm x 18mm	\$267	\$13.35
98.5mm x 20mm	\$283	\$14.15
98.5mm x 25mm	\$317	\$15.85



Flexural Strength Testing per ISO 6872 Specification.  
IPS e.max is a registered trademark of Ivoclar Vivadent.

## BioZX2 Zirconia

BioZX2 is high strength, high translucent zirconia for those cases where strength is the first priority and esthetics are still important. BioZX2 is made in Germany under exacting conditions by Dental Direkt, one of the world leaders in the processing and fabrication of zirconia ([www.dentaldirekt.com](http://www.dentaldirekt.com)).

BioZX2 is perfectly suited for monolithic crowns and bridges of any span, and can be cutback for use with the layering technique.



### Available in 98.5mm discs.

Size	Price	Approximate Unit Cost (based on 20 units per disc)
98.5mm x 10mm	\$142	\$7.10
98.5mm x 12mm	\$167	\$8.35
98.5mm x 14mm	\$183	\$9.15
98.5mm x 16mm	\$200	\$10.00
98.5mm x 18mm	\$217	\$10.85
98.5mm x 20mm	\$233	\$11.65
98.5mm x 25mm	\$267	\$13.35

## DCS - Desktop Containment Unit

**Eliminates Odor, Vapors and Particulate Matter migrating through-out the workplace.**

- Place your system inside the DCS
- Illuminated work area
- Small footprint friendly
- Can operate 24/7 allowing beakers to be left inside contained area
- Designed for quiet operation
- Can be ducted to the outside

**DCS16-GE-M1** includes smaller containment stand (14.5”h x 24”w x 17”d)

**\$1135<sup>25</sup>**  
DCS16-GE-M1

**DCS16-GE-M2** includes larger containment stand (24”h x 29”w x 20”d)

**\$1230<sup>25</sup>**  
DCS16-GE-M2



## JetStream - Small Compact Dust Collector

**Rapidly becoming the standard for small engravers the JetStream offers significant advantages and performance over the competition.**

- Maximum airflow 110 CFM with a maximum pressure of 90” W.C.
- 1 micron pre-filter bag, easily replaced
- 99.97% effective HEPA filter on particles larger than 0.3 micron
- Built in silencer reduces noise to 49 dBA
- Sleek, modern and easily portable
- Full range speed control
- CSA/UL approved motors with circuit breaker and thermal cut-out protection
- Remote on-off communication
- Optional 6 Gallon (26 l) Knock out bucket for larger debris



**\$695<sup>00</sup>**

JS-1200

## FOREDOM<sup>®</sup>



### FOREDOM K.1030 Portable Micro Motor

This all new design has a pocket sized control unit which features a digital speed display at the top so it's easy to read and monitor the speed, which is easy to control with the dial control on the front.

**It operates at a speed of 0 – 30,000 rpm.**

**\$595<sup>00</sup>**

Item: K.1030

## BioTEXT® Wipes

Universal usage disinfectant. May be used on vinyl, plastics, textiles, keyboards, telephones, tablets, smart phones, touch screens, painted surfaces and counters.

- Bioload resistant anti-microbial
- Economical to use
- Multi-purpose
- Non-corrosive
- Non-staining
- User friendly

**\$138<sup>00</sup>**  
 TXW1802

Intro package of (8)  
 & 2 canisters



## Mouth Off™

Anti-microbial for sports mouth-guards, orthodontic appliances, snore guards and dentures. Available in the USA.

- Easy to use spray
- Prevents build up of bacteria and fungi on oral appliances.
- Protects your investment
- Protects you from infections
- Excellent portability and compliance
- Leaves you with a fresh taste in your mouth



## Edentulous Tray Aways Disposable Impression Trays

Edentulous Tray Aways® are single-use, perforated edentulous impression trays that can be easily drilled for implant impressions. These unique and economical trays allow for maximum material retention and accuracy without the use of an adhesive. These high impact styrene trays provide an excellent impression with outstanding details.



**5+1**

keystone  
 Industries

**\$38<sup>95</sup>**  
 0921704  
 Edentulous Assorted pkg (12)

## Mixing Tips \$26<sup>95</sup>

- YM100YB Yellow (50/Pkg.)
- YM110PB Pink (50/Pkg.)
- YM120GB Teal (50/Pkg.)
- YM130PB Purple (50/Pkg.)



## Bibs \$28<sup>95</sup>

- 23-1210 Blue (500/Pkg.)
- 23-1211 Green (500/Pkg.)
- 23-1212 White (500/Pkg.)
- 23-1213 Rose (500/Pkg.)
- 23-1214 Yellow (500/Pkg.)
- 12-1215 Lavender (500/Pkg.)

All Colours!

## Face Masks



Item: #5614  
**\$5<sup>25</sup>**

**5+1**

### NITECH®

#### RONCO NITECH®

- Patented blend of premium quality synthetic materials
- Powder free



### NITRILE

#### RONCO NITRILE

- Diamond textured for superior grip
- Powder free



#### RONCO TOUCH

- Premium quality nitrile in pink colour
- Powder free



### LATEX

#### RONCO SILKTEX

- Superior barrier protection with excellent elasticity
- Available as powder free and lightly powdered



**10+1**

**NEW!**

FROM DETAX THE MANUFACTURER OF  
THE INDUSTRY STANDARD SOFT RELINE  
MATERIAL MOLLOPLAST B

**DETAX**

## Permaliner®

- Shorter polymerization time (1 hour only)
- Lighter natural colour
- Excellent bond, adhesive no longer required

**\$85<sup>95</sup>**

02092 - 30 g


**\$168<sup>95</sup>**

02146 - 70 g



## NU-DENT Denture Cleaner

For the safe and effective removal of coffee, tea, tobacco, tartar and calculus without damaging the appliance.



**64 cleanings per bottle**  
Cleans in less than 25 min.



943230

**\$4<sup>20</sup>**

## Dental Vault

A dental appliance storage container that can be used to clean and store retainers, night guards, snore guards, dentures and partials.



**FREE** Dental Vault  
with a purchase of a  
NuDent Case (48)

Looks like  
a Candle!



05585 Denturevault Dental  
Storage & Cleaning Container (1)

## Duplistar 70

## Duplistar 90

**Ruthinium®  
Group**



2180008

Special

**\$85<sup>00</sup>** Reg.\$115  
1,5 + 1,5 Kg

2180010

Special

**\$176<sup>25</sup>** Reg.\$235  
5 + 5 Kg

2180009

Special

**\$85<sup>00</sup>** Reg.\$115  
1,5 + 1,5 Kg

2180011

Special

**\$176<sup>25</sup>** Reg.\$235  
5 + 5 Kg

A-Silicons for control masks preparation, occlusal registration keys, duplication of full and partial dentures, masks for temporary prosthesis preparation, ceramic moulding, partial counterdies in flask

- excellent details duplication
- dimensional stability (< 0,2%)
- high elastic return
- absence of deformation in flask
- Perfect compatibility with plasters
- Drilling precision

## Sympro

SYMPRO is a compact, high-performance cleaning unit for dentures, orthodontic appliances and splints.



**\$1240<sup>00</sup>**  
6500-0050

**FREE SYMPROfluid**  
Universal with a  
purchase of SYMPRO

6500-0600



## MT3

The MT3 wet trimmer is impressive due to its exceptionally powerful motor. With the MT3 Renfert provides high-performance trimming at a low price.



**\$1345<sup>00</sup>**  
1808-1000

**FREE Klettfix Discs**  
(5 pk) with a purchase  
of MT3 Unit

Klettfix 80/120



## Twister

Twister is a compact mixer with easy-to-use standard functions and powerful motor technology.



**\$2185<sup>00</sup>**  
1826-0000

**FREE Bowl or Stand**  
with a purchase  
of a Twister Unit

Any Size Bowl!



## Waxlectric II

Waxlectric II are the professional instruments in the electric wax knife sector. The application concept concentrates on safety, precision and user-friendliness.



2157-1000

**\$580<sup>00</sup>**

**2 FREE Sculpting Tips**  
with a purchase  
of Waxlectric II



## Waxprofi

Your individual skill is required to give a denture freshness and vitality. This Waxprofi is required to implement this efficiently. Enjoy the pleasure of how easy setting-up has become.



1440-1000

**\$481<sup>40</sup>**

**FREE Wax with**  
a purchase  
of Waxprofi



## Basic ECO

- Compact and convenient sandblasting unit
- Precise working procedures and great cost efficiency
- Innovative Perfect View Concept unique LED technology



2949-4025

**\$995<sup>00</sup>**

**2 FREE Blasting**  
Nozzles with a  
purchase of  
Basic ECO



## Ceramicus

Porcelain brush made from a special, synthetic microhair structure.



**2+1**

6 Brush Sizes:  
(1, 2, 4, 6, 8, 1/0)

## lay:art Style

lay:art style natural hair brushes provide premium quality for the highest demands.



ITEM	DESCRIPTION	REG.	SALE
17250008	Lay:art size 8 slim	\$126.90	<b>\$104.97</b>
17250006	Lay:art size 6 slim	\$93.45	<b>\$77.20</b>
17250004	Lay:art size 4 slim	\$55.95	<b>\$46.25</b>
17250002	Lay:art size 2	\$55.95	<b>\$46.25</b>
17250001	Lay:art opaque brush	\$55.95	<b>\$46.25</b>
17250000	Lay:art brush colour	\$46.20	<b>\$38.15</b>
17251000	Lay:art Style Set 1x #2, 1x Opaque, 1x #1	\$72.70	<b>\$59.95</b>

## Takanishi

Very high quality synthetic bristle brush which has been well proven for many years for the optimal porcelain build-up.



**2+1**

5 Brush Sizes:  
(1, 2, 3, 4, 5)

## GEO Classic modeling waxes

GEO Classic modeling waxes are ideally suited for building up wax technique.

**2+1**



**\$31<sup>50</sup>**

75 g (2.63 oz.) beige-opaque 4920100  
75 g (2.63 oz.) grey-opaque 4920200  
75 g (2.63 oz.) mint-opaque 4920300

75 g (2.63 oz.) beige-opaque 4950100  
75 g (2.63 oz.) grey-opaque 4950200  
75 g (2.63 oz.) mint-opaque 4950300

## GEO Avantgarde modeling waxes

There are special requirements when using the electrical wax up method. These are fulfilled by GEO Avantgarde waxes. The ideal modeling characteristics are maintained and provide precise results.

**2+1**



**\$31<sup>50</sup>**

75 g (2.63 oz.) beige-opaque 4970100  
75 g (2.63 oz.) grey-opaque 4970200  
75 g (2.63 oz.) mint-opaque 4970300  
75 g (2.63 oz.) tooth color-opaque 4970400

75 g (2.63 oz.) beige-transparent 4971100  
75 g (2.63 oz.) blue-transparent 4971200  
75 g (2.63 oz.) mint-transparent 4971300

## Locking Pin Easy-Snap E

Perfect locking pin system for secure grip of the prosthesis. Its small size permits versatile use.

- New mechanism provides secure hold for the prosthesis
- Noticeable snap when opening and closing gives patients more security
- Small size for every situation
- Can be used even in difficult locations
- Simple fitting possible in three variants



Bonding/polymerisation of locking pin Easy-Snap E

The primary construction is always prepared in the same manner for every application. This reduces the need for clarification and accelerates the manufacturing process.



Simple modelling of the primary construction according to the model situation. A set-up for exactly determining the position of the lock is always advantageous.



The hole is measured using the centre drill. Slipping of the Diatit multidrill is thus prevented.



Drilling is simplified using drilling and milling oil. Drill only once, as this prevents the formation of an oval hole.



The hole is closed with the wax.



The wax is removed again by hand with the Rapidy 2.0. In doing so, a slight hollow is created which positions the investment material model exactly.



The model is completely prepared for doubling. If you are working in the withdrawal procedure, omit this step.



Cast the investment material model using the brented doubling system. The high level of design accuracy of the investment material facilitates further processing.



The wax guide is fixed lightly in the cast mould. The framework is modelled with a minimum thickness of 0.5 mm.



The wax guide is removed and an image is created for the locking pin.



The exit hole for the locking pin is easily recognisable on the opposite side.



The locking pin is fixed in the correct position for completion with Qu-resin.



The prosthesis is completed.



Locking Pin Easy-Snap E  
1 piece  
REF 440 0N65 8



Modelling pin E  
1 piece  
REF 440 0065 6



Ceramic spacer E  
2 pieces  
REF 440 0065 7



Device for waxing-on  
1 piece  
REF 440 0066 1



## Abraso-Sil Acryl

This buff consists of a nonwoven fibre fabric in the centre between two layers of silicone-coated cotton fabric.

On the outside there are two rows of bleached Chungkink bristles. This brush absorbs polishing paste or pumice and releases it slowly for efficient prepolishing.



Ø 80 mm  
1 piece

**8 + 1 FREE!**

## ROPAK/Kompaktopaker UV



**Ready to use Light Cure Opaker!**

The ready-to-use alternative for convenient coating of CoCr objects.




Apply Ropak Kompaktopaker with the integrated brush directly onto the clean metal surface.

Reg. \$139.29 Ropak Kompaktopaker (Pink) 20ML 54000133  
**\$99<sup>50</sup>**

Reg. \$76.57 Kompaktopaker Zahnfarben (Tooth Colour) 10ML 54000105  
**\$56<sup>00</sup>**

## Breident Diagen Turbo

Diagen turbos are very aggressive sintered diamonds that work effectively on zirconium oxides such as Procera™, Emax™ and all other ceramic materials. Each stone is molded to provide optimum grinding properties and cuts effectively down to the last layer.



- |   |   |   |
|---|---|---|
| <br><b>\$85<sup>00</sup></b> | <br><b>\$60<sup>00</sup></b> | <br><b>\$65<sup>00</sup></b> |
| Cylinder<br>Ø 6.5 x 13 mm<br>2 pieces<br>REF 340 0017 0   | Cone<br>Ø 3.5 x 11 mm<br>2 pieces<br>REF 340 0015 0   | Inverted cone with recess - Ø 12 x 6 mm<br>1 piece<br>REF 340 0024 0  |
| <br><b>\$64<sup>00</sup></b> | <br><b>\$64<sup>00</sup></b> | <br><b>\$85<sup>00</sup></b> |
| Lens<br>Ø 22 x 2 mm<br>1 piece<br>REF 340 0021 0  | Disc<br>Ø 22 x 2 mm<br>1 piece<br>REF 340 0022 0  | Disc<br>Ø 22 x 4,5 mm,<br>1 piece<br>REF 340 0019 0   |

**6 + 1 FREE!**  
*OF EQUAL OR LESSER VALUE!*

## Universal Screwdriver Set

Screwdriver set for 98% of all screws available on the market. To be inserted into the torque ratchet, adjustable from 10 to 40 Ncm. This way screws can be turned in correctly and safely.

Universal screwdriver set with instruments  
REF 310 0001 2

**\$1995<sup>00</sup>**  
Reg. \$2233.00



## Torque Ratchet

Torque adjustable from 10 to 40 Ncm.

**\$350<sup>00</sup>**  
Reg. \$508.00



Torque ratchet  
REF 330 0115 5

# Ruthinium CAD/CAM DISCS




PMMA discs for fabrication of longer lasting crown and bridges, manufactured with high molecular cross linked acrylic which guarantees high detailed precision and resistance.

## Advantages:

- High milling reproduction details
- Absence of cracks even on thinner layers
- Resistance to abrasion
- Fluorescence
- Burns out completely
- Chromatic stability and absence of porosity
- Biocompatibility



	RUTHINIUM DISC SINGLE LAYER	RUTHINIUM DISC MULTILAYER
<b>Shades</b>	A1 - A2 - A3 - A3.5 - B1 - B2 - C3 PINK - TRANSPARENT - RX	A1 - A2 - A3 - A3.5 - B1 - B2 - C3
<b>Code</b>	7900087	7900088
<b>Diameter</b>	98.5MM	
<b>Thickness</b>	16 mm - 20 mm - 25 mm	
<b>Price</b>	<b>\$85.00</b>	<b>\$125.00</b>



**Implant Overdentures NEED TO BE REINFORCED!**

**HOT DEAL**

Buy the FiberForce kit for \$538 and get  
1 pack of FiberForce UD and 1 pack of  
FiberForce BRAID 1mm for FREE!

### FiBER FORCE® starter kit includes:

- 2 - packs of FiBER FORCE mesh (pink)
- 1 - pack of FiBER FORCE 2mm braided rope (pink)
- 1 - syringe of CST Flow pink light cure resin (3ml)
- 100 - sheets of pink spacer wax (0.5mm)
- 100 - patient brochures (with cardboard stand)
- 1 - EZ VAC vacuum forming unit



**"ALL ON CST®"**

**HOT DEAL**

Buy the CST starter kit and get  
1 Refill of cst for FREE!

### CST® complete kit includes:

- 2 - 450mm CST® Pink Hybrid Compressible Fiber (1:6)
- 2 - 300mm CST® Pink Hybrid Compressible Fiber (1:4)
- 1 - 3ml syringe of CST® Flow Pink Resin with tips
- 1 - 5ml bottle of CST® Bond Resin
- 4 - CST® Temporary Support Pillars



## Denture Brushes



Purchase a box of 50 and receive 12 FREE Denture Baths!

## Acry-Tray



0800001

Buy 1 package & receive FREE 1Lb. Alginate with Spatula!

FREE 1 lb. Alginate #2060001 & 1 Spatula #904SA

\$69<sup>95</sup>

## Elite HD

4 + 1 FREE!



## Zhermack Alginates



5 + 1 FREE!

## Tokuso Rebase II Kit

FREE BUR!  
\$24<sup>value</sup>  
#5410.060



Buy 1 package & receive a FREE Carbide Bur!



20506 Tokuso Rebase Kit  
\$169<sup>50</sup>

## Diamond D - Heat Cure Kit

Buy 5Lb. & receive 25 FREE Denture Baths!



1013065  
5 Lb. Diamond D H/C  
Original P + L

## Easy Soft Economy Kit

FREE Trimming Bur!  
#7275.060



ES0502 - Easy Soft Economy Kit

## Denture Bath

10 + 1 FREE!



## Molloplast B 170g.

FREE Trimming Bur!  
#7275.060



03001



**8<sup>th</sup> ANNUAL SPRING CLASSIC GOLF TOURNAMENT**

**The Briars Golf Club**  
127 Hedge Road, Sutton West, ON L0E 1R0  
Friday, June 3rd, 2016 @ 1pm Tee-off

The Briars Golf Club is located just one hour from Toronto on the southern shores of Lake Simcoe. As you wind your way along historic Hedge Road, you immediately feel relaxed and understand the allure of this location. Two Canadian design legends have combined to bring you a delightful dose of golf at its best. Stanley Thompson carved out the original nine holes in 1922. Robbie Robinson integrated a second nine in 1971, creating a seamless routing that may look fairly innocuous on paper, but will be appreciated for its subtle charms. The Briars is the epitome of the old-style course.

**To Register Call: 1.800.268.4442**

**\$140/person**  
\*Includes 18 Holes,  
Cart, Dinner, Prizes & More!



# WHY PAY MORE?



M82-C-Single Workbench



LAB-21 E2 | LAB-21 D | LAB-21 GC

Plaster Bench (3-Unit Combo)



M-93 Double Workbench



M-SAM

## Good Quality... Great Value!