

3D SYSTEMS

Introducing...

NEXTDENT 5100
Redefining Digital Dentistry

HIGH-SPEED DLP PRINTER


Revolutionary high-speed figure4™ 3D printing technology combined with the industry's broadest portfolio of dental materials addressing multiple applications, resulting in unparalleled accuracy, repeatability, productivity, and total cost of operation.



Nextdent-5100

**Shown with optional Pedestal*

For more information contact us at 1-800-268-4442 or visit website: www.CentralDentalLtd.com




Cadstar Neo and Neo Pro

Structured White Light Lab Scanners

Precision-Made in Austria

Neo - Starting from \$14,750
Neo Pro - Starting from \$15,750



CS.NEO
CS.NEO^{PRO}

Software Features

- Smart Integration - Manual or automatic rescan function to gain even more scan data from hard to capture areas.
- Guided Workflow function - for hassle-free scanning and expert mode with more options and customizable scan strategies for advanced users.
- Click-less workflow - Click-less workflow with full automatic alignment and automated work steps.
- Seamless workflow - integration with Exocad
- No Annual License Fees




Smart operating

3D Mouse and Touch Screen
Ready Navigation.



High resolution

HR scan data with sharper edges
and cleaner details.



Reliable results

Provided by a high precise glass calibration body and a calibration control system. Environmental sensors are giving feedback to the user when it's time for a recalibration.



A tool to increase productivity

Speed up your daily work by reducing scanning steps with the quarter tray adapter.



White LED light source

To obtain as much scan data as possible, CS.Neo uses the full spectrum of light.



TRU color

Hand-drawn color markings on the model are neatly and sharply captured, and prepared in Exocad for the design process. This function is especially helpful for the digital manufacturing of partial frameworks.

NEXTDENT BIOCOMPATIBLE MATERIALS

THE NEXTDENT 5100 ADDRESSES MULTIPLE APPLICATIONS WITH A BROAD RANGE OF NEXTDENT MATERIALS

<p>Pending Health Canada Approval</p>  <p>Denture 3D+ 3D print resin for the manufacturing of removable denture bases</p>	<p>Pending Health Canada Approval</p>  <p>C&B MFH Micro filled hybrid 3D print resin for the manufacturing of long-term temporaries</p>	 <p>Try-In 3D print resin for the manufacturing of try-in devices</p>	 <p>Tray 3D print resin for the manufacturing of individual impression trays</p>	 <p>SG (Surgical Guide) 3D print resin for the manufacturing of dental surgical guides</p>
 <p>Model 2.0 3D print resin for the manufacturing of prosthodontic and orthodontic models</p>	 <p>Gingiva Mask 3D print resin for the manufacturing of gingiva masks on dental models</p>	 <p>Ortho Rigid 3D print resin for the manufacturing of dental splints</p>	 <p>Ortho IBT 3D print resin for the manufacturing of Indirect bonding trays.</p>	 <p>Cast 3D print resin for the manufacturing of castable parts</p>

VHF milling and grinding mills - "Dental from Deutschland"

VHF Camfacture AG has been making milling and grinding machines under private-label for companies around the globe and under it's own brand for thirty years. VHF dental milling and grinding machines are robust workhorses, designed to reliably withstand the stresses created during milling and grinding; especially of harder materials like chrome-cobalt, glass ceramics and PMMA.

VHF offers a fully integrated solution with validated work-flows. Companies in the VHF group manufacture all components of the machines including the internal electronics, the milling and grinding tools, and the dental CAM that define the milling and grinding strategies. This ensures that outputs from the machines are exactly as expected and desired, making manufacturing processes predictable.

THE NEW K5+ IS A 5-AXIS MACHINE DESIGNED FOR PRECISION AND LONGEVITY



- The K5+ b-axis will rotate up to 35 degrees, making the milling of thicker elements like denture base a snap. The higher rotation range allow for milling of discs up to 40mm thick.
- Powerful and robust- mills the toughest materials, including CoCr.
- New direct disc technology allows for one-handed clamping of disc.
- vhf DentalCam7(included) optimizes CAM strategies for the fastest possible machining.

Lease From \$ 928.00/mth
5 YEARS TERM

N4 - WET GRINDING OF GLASS CERAMICS AND TITANIUM

The N4 is a compact 4-axis wet grinding machine suited for the production of glass ceramics, zirconia and titanium. The N4 has a fully integrated liquid cooling system, simplifying the overall equipment needs and setup.

Lease From \$742.43/mth
5 YEARS TERM



Diamond D Denture Base

- Handling and finishing properties are second to none
- High tooth to base adhesion
- Dead on shade matches to industry standards
- Exceptional flexural strength and impact resistance
- Smooth packing consistency
- Better fit

Buy 25lbs of Diamond D denture base resin Get 4qt. of liquid FREE!



Buy 5lbs of Diamond D denture base resin Get 1qt. of liquid FREE!

Tecnoflask Kit

- Does not require a specific resin
- Increased removal of residual monomer
- Reduction in energy consumption
- Cleaner and faster than the traditional
- High quality denture



Buy a Tecnoflask Kit Get 1 qt. of Ultra Sep FREE!
#1009200

keystone industries

ENAMELITE Low-Fusing Ceramic Spray Glaze Glaze ceramics like never before!

- Reduces labor cost, increases efficiency
- Spray a thin, uniform, & accurate coat IN SECONDS
- Faster drying than brush-on applications
- Compatible with ALL dental ceramics



Buy 2 Receive 1 FREE

#5100179

KEYPRINT print resins – “Juice for Your Printer”

Need a print resin for your open printer? Keyprint resins from Keystone Industries will do the trick.

- Keymodel beige - \$325 per kg.
- Keyortho model grey - \$325 per kg.
- Keymask pink : soft tissue model resin - \$285 per kg.
- Keycast: casting resin - \$535 per kg.
- Keyguide: Surgical guide resin - \$500 per kg.

Now available for sale in Canada



Fiber Force Reinforcement & CST frameworks

- Replace cast and milled bars in fixed hybrid implant supported dentures
- Innovative new fiber design and simple fabrication procedure = fiber framework fabrication in 30 minutes

Q2 SPECIALS

- Buy Fiberforce starter kit and get 1Pkg. UD & 1Pkg. Imm Rope Free
- Buy 6 Pkg. mesh get 1Pkg. Free
- Buy 6 Pkg any Rope/UD (mix & match) and receive 1 Pkg. of any Rope or UD Free



BRUSHLESS MICROMOTOR NEW!

From Freedom

High torque at all speeds
with a maximum of 50,000 rpm

2 year warranty

Intro price
\$750.00



A7 Plus Articular

Articulator A7 Plus is made of aluminum,
assuring greater stability.
Design that provides greater visibility.

- Semi-adjustable / Arcon Type
- Fixed Intercondylar distance at the average 110 mm.
- Curved condylar guide
- Adjustable condylar guide angle
- Adjustable Bennet Angle
- Central Lock
- Stabilization system of Condylar Guide
- movement with silicone connection
- Pin to support the Upper Frame at the open position



Includes Articulator,
Face Bow & Carrying Case

Face Bow



Lower frame designed to
allow inclination at $\sim 45^\circ$

Bennet and Condylar
Guide angle adjustable


Curved Condylar
Guide Central Lock



BRENT Easy Snap E Locking Pin

- The locking pin Easy-Snap E is used for removable dentures, free-end saddles, partial dentures, telescopic and conical crowns or implant supported dentures.
- The snap rings supported guidance results in a soft snap of the pin during locking in the open or closed position.

ITEM #4400N658
\$89.39



DTK ADHESIVE KIT

- A dual hardening composite adhesive for the fixation of dental attachments

#54001185
\$166.00

Bredent Diagen Turbos - Fine Grit!

Diagen turbos are very aggressive sintered diamonds that work effectively on zirconium oxides such as Procera™, Inceram™ and all other ceramic materials. Each stone is moulded to provide optimum grinding properties and cuts effectively down to the last layer. Both the porcelain bonding material and the synthetic diamonds are carefully proportioned in order to ensure maximum cutting efficiency and minimal heat so that chipping is virtually eliminated.



Bredent Diagen Turbos: Unique Shapes and Coarser Grits

When Bredent learned that technicians were sometimes modifying the shapes of the Diagen Turbo to accommodate their specic needs, they responded with NEW shapes - a pointed cylinder which is perfect for interdental areas and a rounded cylinder which is ideal for nishing collars and chamfers. They also have extended the product range to include a Diagen Turbo with a Coarser Diamond Grit (yellow) which is 20% more aggressive than the regular Diagen Turbos.



Buy 6, get
1 FREE

red & white with crea.lign®

Red & white with crea.lign® is the solution for:

- Implant-supported, fixed or removable permanent restorations
- Natural customization of full dentures
- All substructure materials such as BioHPP® and PEEK
- Complete coating even of full-arch restorations (e.g. Prettau bridges)
- Can also be used with ceramic and non-precious metals

SPECIAL

\$595.00

Reg. \$675.00

CLIGNSETC

bredent group

bre.lux PowerUnit 2

\$3395.00
GREAT VALUE

EASE OF USE

- Easy and fast access to programs
- Programmable to satisfy individual requirements
- Spacious drawer
- Extremely silent unit
- Compact unit

Reliable polymerization

The NEW Innovative LED Technology of Brelux 2 covers the entire relevant spectrum of light, allowing for a more consistent, faster cure.



#14001000

Unsurpassed Polymerization Speed

- 72 lights + full range + rotary plate + Light-Tray

Long service life of the LED lights

- LED Lights rated for 20,000 operating hours or 12 years.



CANDULOR

KMG®

High-gloss Polishing Paste

Universal high-gloss polishing paste for acrylics and metals. It is easily washed off and particularly economical to use

BENEFITS

- Simple to apply
- Suitable with any buffing wheels
- Washes off easily
- Economical to use
- Slightly minty taste



Back by Popular Demand

\$ 17.61
#693113
500ml

ISO-K®

Separating Liquid

Alginate-based separating liquid with a proven success record. It is ideally suited for isolating plaster or stone against denture base materials, i.e. PMMA acrylics of all kinds, both heat curing and cold curing



BENEFITS

- Excellent separation
- Readily absorbed into the plaster
- Forms very fine separating film

\$ 30.96
#559537
500ml

EDENTA

SuperMax® Diamond Trimmers

Sintered diamond with organic bonder for grinding of ceramics and metals. Enables rough trimming without heat build-up in the restoration. Suitable for rough trimming of Zirconia, porcelain, all-ceramics, and pressable porcelain sprues.



SPECIAL Buy 5, GET ONE FREE!

9001-9003 \$28.70
9004-9007 \$35.15

StarGloss® Diamond Polishing System

Comprehensive polishing system specially designed for gently polishing porcelain and all-ceramic restorations (including Zirconia). Natural diamond grit with a condensed, long-mesh synthetic bonder, produce an outstaying surface texture Three-step system eliminates polishing pastes or extra glaze rings.



Blue - Coarse (Pre-Grinding/Bulk Reduction)
Pink - Medium (Smoothing/Pre-Glaze)
Grey - Superne (High-Shine/Post-Glaze)
Wheel: \$38.55
Knife Edge: \$51.05
Flame: \$12.65

SPECIAL Buy 5, GET ONE FREE!

PASCAL

Hydent

Denture Indicator Paste

- Fast - Apply to denture with the touch of a button
- Accurate - Pressure points are clearly defined
- Even application with aerosol spray
- Easily cleaned from denture
- No messy brushes
- Mint flavor
- Contains no CFCs



\$ 31.95
 Plus Free Bottle
 Nudent Denture Cleaner

#05100 - 30g

TALLADIUM

TuffRock

Type IV Die Stone

A synthetic resin die stone formulated with a low expansion of .05% in 1 hr. & up to 5 days & longer. Tu Rock has an excellent hardness and strength without the brittleness (20,450psi).

Colours available:

Blue, Green, Ivory, Peach, Pink, Topaz & White



\$75.95
 (10kg)
 5 + 1 FREE

SHADEWAND

LED SHADE MATCHING LIGHT

Shade matching perfection!

The shadewand is a held, rechargeable LED shade matching light that brings the perfect light to the patients mouth. This 9" shade wand eliminates ambient light variables such as wall colour, clothing bids etc. These light variables effect shade selection. The shade wand uses 5500 kelvin led which generates the perfect light for taking shades.

- rechargeable
- 3 light settings
- 5500 kelvin light
- easy to use

Special
 \$ 169.95
 SAVE 40.00



161

RUTHINIUM

ACRY- PLUS Premium Quality

- 4 Layer Acrylic Teeth for removable dentures.
- Unique Enamel-Dentine cutting has been conceived for the upper anteriors, thus increasing the "Translucent Effect" making the Acry-Plus Teeth Naturally Characterized and Life-Like!

NEW OFFER

PURCHASE \$500.00 OF ACRY PLUS TEETH
 RECEIVE \$75.00 GIFT CARD

PURCHASE \$1500.00 OF ACRY PLUS TEETH
 RECEIVE \$250.00 GIFT CARD

(PURCHASE HAS TO BE ON A SINGLE INVOICE TO APPLY)



(1 x 6)
 \$25⁵⁰



(1 x 8)
 \$19⁸⁰

Twister Mixer



Bonus
Your Choice of
FREE STAND
or Mixing Bowl

Special
\$ 2215.00
Reg. \$2522.65
1826-1000

Waxprofi II



Bonus
1 - 5lb of C/D #2
Baseplate Wax

Special
\$ 435.00
Reg. \$476.40
1440-1000

MT3 Model Trimmer

Bonus
Klettfix 1 pkg.(5)
Trimming Discs
\$ 115.00 (value)



Special
\$ 1299.00
Reg. \$1482.50
1808-1000

Waxlectric II

Bonus
2 FREE Tips!



Special
\$ 440.00
Reg. \$501.00
2157-1000

Sympro Cleaner

Bonus
FREE 2 x 2 L
Sympro Universal
Fluid \$98 (value)



Special
\$ 1285.00
Reg. \$1460.00
6500-0050

Basic ECO

Bonus
5 Kilo Cobra 50
Micron Aluminum
Oxide \$55 (value)



Special
\$ 1040.00
Reg. \$1189.35
2949-4025

STOCK UP FOR THE SPRING • Everyday Best Sellers

Molloplast B
170g.



#03001

Reg. \$281.99
Special \$256.99

SAVE \$25

Lucitone 199
25lb.



SAVE \$100

Reg. 1102.60
Special \$999.00

Edenta



5+1 Burs

C/D #2 Baseplate Wax



SAVE \$5

Reg. \$49.50
Special \$44.50

#0252

Elite HD
(variety)



4+1 FREE

Zhermack Alginate



10 + 2 FREE

Micro-Torch



#04020

Excelent VALUE
\$31.95

Denture Baths



Imprinting Now Available!

BUY 100 Receive 10 FREE

Zetalabor



3 + 1 Available in 2.6kg
5kg & 10kg

EXA-Technique



6 + 1 FREE

Cotton Buffs

4 x 42



#1180170

**Purchase 1 Dozen &
Receive 2 FREE Buffs!**

B12 Brushes



#B12

**Purchase 1 Dozen &
Receive 2 FREE Brushes!**

C/D Repair Liquid



SAVE \$10

Reg. \$72.60 #00323
\$62.60

Al-Cote



Reg. \$127.85 #652502
\$115.06

DENSEP Separator



SAVE \$10

Reg. \$64.80 #555-8
\$54.80

Denture Brushes



#00002

Reg. \$83.99
\$73.99

50 Pkg.

SAVE \$10

11th. Annual Spring Classic

The Briars Golf Club is located just one hour from Toronto on the southern shores of Lake Simcoe. As you wind your way along historic Hedge Road, you immediately feel relaxed and understand the allure of this location. Two Canadian design legends have combined to bring you a delightful dose of golf at its best. Stanley Thompson carved out the original nine holes in 1922. Robbie Robinson integrated a second nine in 1971, creating a seamless routing that may look fairly innocuous on paper, but will be appreciated for its subtle charms. The Briars is the epitome of the old-style course.

\$150/person
*Includes 18 Holes,
Cart, Dinner, Prizes & More!



**The Briars Golf Club
GOLF TOURNAMENT**
127 Hedge Road, Sutton West, ON L0E 1R0
11th. Annual Spring Classic
Friday June 7th. 2019 @ 1.00pm shotgun start.
To Register Call: 1.800.268.4442

EVENTS CALENDAR



EVENTS:

- SDA - March 2nd
- DTA - March 23rd
- Technorama - April 5-6th
- Spectrum Day Saskatoon - April 13th
- Spectrum Day Edmonton - April 26-27th
- Spectrum Day Montreal - October 4th
- Spectrum Day Toronto - October 25th
- Spectrum Day Vancouver - November 9th
- DGO - Annual general meeting and trade show April 26th. Toronto
- Denturist Association of British Columbia May 29th-June 1st. Penticton
- Denturist Association of Alberta September 11-14 Edmonton
- Denturist Association of Ontario- PYP- TBD
- Symposium Quebecois de la Denturologie September 11-14 Trois Rivieres