

NFC⁺

3 tooth lines that will take your breath away. **NFC⁺ by Candulor.**



ZTM Tanja Suppiger

- + Phenomenal Esthetics
- + Top Abrasion Values
- + New Material Formula



Bonartic[®] II NFC⁺



Condyloform[®] II NFC⁺



PhysioStar[®] NFC⁺

www.candulor.com

Swiss Denture Concept

- A Special Training Course For Denturists.

For more info go to: www.swissdentureconcept.com

Check the Events Page on our Website:

www.CentralDentalLtd.com

For the Swiss Denture Course in your area.

GEO Wax Wire

Wax wire on spools for model casting techniques

\$16.95 each

**BUY 2
GET 1
FREE***



NEW!

RENFERT-Scanspray

High end surface matting

- The container size of 200 ml provides an optimal price-performance ratio.
- Maximum precision thanks to very thin layer application and uniform grain size.
- Economical to use due to its high covering power.
- The special valve technology ensures spraying without clogging.



#1731-1000

Reg. \$79.95

\$68.95 each

K-Diamonds

edenta | e

All-ceramic ZrO₂ preparation with diamond instruments using water cooling

A prerequisite for highly durable restorations is material-friendly preparation of the sintered ceramic in order to avoid microcracks and ceramic splitting off. Large surface areas should no longer be prepared, but only essential, minor fitting adjustments using special K-Diamond rotary instruments with water cooling. The 3-phase multilayer technology in combination with a newly developed solid nickel matrix guarantees a high abrasive capacity with maximum service life. The grit size of the diamond instruments has been specifically designed to avoid removal of ceramic particles to ensure that there is no risk to the long-term durability of the allceramic restoration.

K-Diamonds

200.000
5



Fitting the sintered zircon oxide framework

By trimming the high spots with diamond instruments.
For trimming and adjusting sintered /hot isostatically pressed (HIP) zircon oxide using a laboratory turbine and cooling with water.



Shank	L mm	ISO	Order No.	5,5	8,0	11,5	10,0	8,0
FG		•• K806 314 263 514...	KF369.314...	025				
		•• K806 314 263 504...	KC369.314...	025				
		•• K806 314 263 494...	KUF369.314...	025				
		•• K806 314 697 514...	KF801L.314...		014			
		•• K806 314 697 504...	KC801L.314...		014			
		•• K806 314 697 494...	KUF801L.314...		014			
		•• K806 314 198 514...	KF856.314...			016		
		•• K806 314 198 504...	KC856.314...			016		
		•• K806 314 198 494...	KUF856.314...			016		
		•• K806 314 167 514...	KF859L.314...				010	
		•• K806 314 167 504...	KC859L.314...				010	
		•• K806 314 167 494...	KUF859L.314...				010	
		•• K806 314 290 514...	KF879.314...					014
		•• K806 314 290 504...	KC879.314...					014
		•• K806 314 290 494...	KUF879.314...					014
		•• K806 314 141 514...	KF881.314...					016
		•• K806 314 141 504...	KC881.314...					016
		•• K806 314 141 494...	KUF881.314...					016

\$37.95
5 / Pkg

HARVEST DENTAL

ZCAD Wax Press
Milling wax disc for
pressable and lost
wax casting.



Item #1006

\$29.95

PRECISION:
OPTIMIZED FOR INTRICATE
HIGH SPEED MILLING.


Benefits

- Stable and strong.
- Machine to total marginal accuracy.
- Burns-out clean and residue free.
- Work with water and dry Cam systems.

With a disc size of 98.5 x 18mm, they are compatible with all open CAD/ CAM systems that accept milling discs with a diameter of 98 mm. Includes step.

**ZCAD PolyCast - Polymer
Milling Disc**

Fiber-free acrylic polymer
for milling accurate
cast models.



Item #1010

STABILITY:
HIGH RESISTANCE AND CLEAN
BURNOUT.

\$52.95

Benefits

- High resistance and dimensionally stable.
- Allows optimal thickness of walls.
- Burns out clean and residue free.

The crystal clear, high translucency of the disc makes it ideal for milling surgical guides, offering excellent visibility for easy implant placement.

TUTTNAUER

EZ9P & EZ10P Autoclave with Printer



Special
EZ9P **\$5185.00**

Special
EZ10P **\$5795.00**

PRESTIGE

Prestige Classic 2100
Portable Autoclave



Cost Effective
Solution for
Sterilization

\$1250

#210900

ER METALS Lab Workbenches

- Metal surfaces are corrosion resistant
- Painted by electrostatic powder method with many colour options
- Many options for material of work-top
- Retractable air nozzle and easy to use gas tap
- Vacuums with dust bags are an option
- Arm support and chromium work-sheet for the worktop
- Easily accessible on/off switch with electrical outlets
- Two-styles of light fixture available
- Single-station, single-set of drawers and vacuum

Starting from
\$2,615



#M82



#M93



#M96



#M82a

GARRECO

CAD-scan
Scannable Die Stone - Extra-Hard
Set Time: 11-12min.
Comp Strength: 18,500psi
Expansion: 0.06-0.09%



Regular \$101
 #1462550
 Green (50lbs)

Special
\$95.95

RUTHINIUM

AcryRock Acrylic Teeth


Available in ALL 16 Vita Shades
 Large Mould Range
 - 32 Upper & 19 Lower Anteriors

**THE BEST VALUE
 IN THE MARKET
 FOR ECONOMY
 TEETH!!!**



\$3.95

/card



DuraFlex Partial



DuraAcetal Framework

myerson FlexPress™

automated digital injection system

Reliable, Repeatable Results

\$100 FREE MATERIAL

FlexPress Professional Kit

Reg. \$5503.65

Special Price

\$4995.00

Plus \$100 Free Material

#FP-0101

FlexPress Advanced Kit

Reg. \$4702.95

Special Price

\$4195.00

#FP-0100

FlexPress Machine Only

Reg. \$4421.76

Special Price

\$3995.00

#FP-0110



Key features

- Integrated heating unit both melts thermoplastic material and preheats injection flask.
- Digital controller ensures highly accurate temperatures.
- Melt Temperature, Melt Time, and Hold Time for myerson DurAcetal[®], myerson DuraFlex™, and myerson VisiClear[®] are pre-programmed.
- Repeat an injection profile with just the touch of a button.
- Program and save your own settings (3 maximum)
- USB port allows easy upgrade of software
- 3-year warranty
- Made in the USA

*Requires hookup to air pressure, minimum 90psi (6.2 bar)

myerson
the tooth company
more to myerson than teeth

myerson HotShotElite™

handheld clasp injector

make clear clasps in under 30 minutes



Professional

Reg. \$2494.95

Special Price

\$2195.00

#EG-0102

Advanced Kit

Reg. \$2111.00

Special Price

\$1895.00

#EG-0101

Gun Only

Reg. \$1726.95

Special Price

\$1550.00

#EG-0103

Basic

Reg. \$1918.95

Special Price

\$1695.00

#EG-0100



GC FUJIROCK® EP

Synthetic Die & Model Material

GC FUJIROCK EP is the ultimate stone for producing remarkably fine/smooth dies and working and study models. The superior quality of this Type IV Dental Die Stone is ideal for all kinds of prosthetics with high precision, outstanding edge hardness and high pressure stability. This product features a longer pour time and a faster set time. It resists chipping, exhibits high compressive strength, lowers expansion and yields a hard, glossy surface. One of the world's most popular dental die stones, GC FUJIROCK EP, is available in seven colors, which meets the requirement of all demands for aesthetic and high quality dental die stones.



**BUY 3
GET 1
FREE***

NATURE-CRYL® SUPER HI IMPACT

Ultra High-Impact, Denture Base Resin

- Ultra high impact to resist breakage or chipping
- Easily finished and polished
- Superior physical properties
- Accurate fit
- Convenient 2-hr heat cure
- Overnight cure
- Available in 7 different shades



NEW!

**Purchase 5lb
Powder and receive
1 Qt Liquid FREE!**

Intro Special

\$147.50
5lb Powder



#43106
Stone Chuck
1/4 (right side)



\$4.95
Special

#42150 - Magnetic Bur Holder



\$15.75
Special

#1150080 - Mixing Bowl



\$2.95
Special

#26207- Wax Cup

Reg. \$9.95
\$7.95
Special




#26101
Bunsen Burner

Reg. \$33.95
\$25.95
Special



#D11A
Silicone Mixing Bowl
With Lid

Reg. \$18.95
\$9.95
Special



#41301
Boley Gauge

\$19.95
Special

\$7.75
Special



#43114
Taper Chuck Steel
Left Side
#43114 - Right Side

The *Miracle* of New Ultrathin Nitrile!

NITRILE POWDER-FREE EXAM GLOVES

A truly superior performance you would never expect from a nitrile glove!
 Optimal stretchability, just like latex! Long-lasting comfort, durability and decreased hand fatigue. Beaded cuff snugs wrist, eliminates roll-up, minimizes entry of liquids. Thickness: 2.5 mil. Finger-textured finish for precision gripping. Latex-free, non-allergenic.

Color: Teal Green.



200
Gloves / Box

Double the
Normal Amount! (100/box).

Intro Special
\$16.95 / Box

8+2

CODE/SIZE
92395-XS
92396-S
92397-M
92398-L
92399-XL

MICRO-KLEEN 3

Is ready to use, broad spectrum surface & equipment disinfectant cleaner. Micro-Kleen3's fast "Accelerated Kill Time Formula™", will effectively clean and disinfect surfaces such as: Dental equipment, operating lights-tables, impression materials prosthetics, counter tops, sinks, reception areas and all "HIGH CONTACT, HIGH TOUCH" environmental health care surfaces.



(#MKS12-024-1) 24oz. Sprayer bottle **\$9.50**

(#MKW06-100) Case of 6-Wipes (100)/Roll **\$77.70**

(#MKP04-128-1) 3.8L bottle **\$24.95**

INVENTORY CLEARANCE SALE!

- #2060001 Alginate Dust Free(300 only) Regular- \$9.95 **Special \$6.00/Bag(500gm)**
- #42150 Magnetic Bur Holder(25 only) Regular- \$18.95 **Special \$14.00/ea.**
- #C0335 C/D Economy Baseplate Wax 5LB. Box(200 only) Regular- \$29.95 **Special \$22.95**
- #AF400 DM Quatro Air Purifier(3 only) Regular- \$1202.00 **Special \$895.00**
- #62-II Handler Super Sucker Dust Collector(2 only) Regular- \$719.80 **Special \$550.00**
- #903090 Success System complete(1 only) Regular- \$9046.07 **Special \$5500.00**
- #42060 BDM Model Former Small(29 only) Regular- \$6.15 **Special \$2.00/ea.**
- #42063 BDM Model Former Medium Lower(10 only) Regular- \$6.15 **Special \$2.00/ea.**
- #42061 BDM Model Former Small Lower(only 71) Regular- \$6.15 **Special \$2.00/ea.**
- #78270 BDM #61 Plate Shear(30 only) Regular \$24.80 **Special \$15.00**
- #F-201 Corident Twin Pins(41) Regular \$99.95 **Special \$35.00/pkg(1000)**
- #F-200 Corident Plus Pin Dowel Pins(37) Regular-\$69.95 **Special \$35.00**
- #03375 Detax Easyform LC Paste(13 only) Regular- \$106.75 **Special \$69.95**
- #SW725-50 Master Bond Brown rubber points (20 only) Regular- \$29.95 **Special \$9.95/box**
- #MTD10C Model Trimmer Disc 10inch (3 only) Regular- \$650.00 **Special \$300.00**
- #MTD12C Model Trimmer Disc 12inch(5 only) Regular- \$825.00 **Special \$350.00**
- #IV-811 Quatro I-Vac dust collector (1 only) Regular- \$1395.00 **Special \$1095.00**
- #305-311 PTC Matrix waxing instrument (15 only) Regular- \$83.95 **Special 49.95**
- #7009011- Preat perma mesh sheets(11 only) Regular- \$142.70 **Special \$65.00**
- #200TRD-A Pro Retainer box(pkg 12) assorted (158 pkg 12 only) Reg- \$7.20 **Special \$3.60**
- #1670140 Hatho cotton brush with wood center(15 only) Reg- \$14.45 **Special \$8.00**
- #7000371 Work Pan Blue (1500 only) Reg- \$2.20 **Special \$1.50 ea.**
- #CSU3 Tuttnauer ultrasonic cleaner 3 gallon (2 only) Regular- \$1285.00 **Special \$895.00**
- #CSU3H Tuttnauer ultrasonic cleaner 3 gallon with heater(2 only) Reg-\$1450.00 **Special \$995.00**
- #CS250 Digital Comunicator II(2 only) Regular- \$650.00 **Special \$375.00**
- #ENZOU – Enzo Universal Articulator(9 only) Regular- \$275.00 **Special \$175.00**
- #08842- Whip-Mix Vac-U-Mixer 1200ml.(20 only) Regular- \$227.60 **Special \$85.00**
- #31315- Whip-Mix Lean Rock 50lb Ivory(5 only) Regular- \$108.50 **Special \$45.00**
- #2524-6000 Renfert vortex Compact dust collector(1 only) Regular- \$953.40 **Special \$650.00**
- #2954-1005 Renfert Quatro 2 tank Sandblaster(1 only) Regular- \$1797.97 **Special \$950.00**
- #2921-1050 Renfert Silent Dust Collector-(1 only) Regular- 2152.00 **Special \$750.00**
- #550 Handler Nu Style Port-A-Vac (2 only) Regular- \$390.40 **Special \$275.00**