



# Refresh.

verb  
re-fresh \ri-'fresh\

: To give new strength or energy to;  
To reinvigorate the Customer experience  
By exceeding expectations  
At every level, all the time.



Central Dental. Now a proud supplier partner of Harvest Brands.

**C/D BITE RIMS**



**\$75<sup>60</sup>** 47415  
4lbs.

**C/D REPAIR LIQUID**



**\$64<sup>00</sup>**

00323  
1 Quart

**C/D DENSEP**



Special  
**\$49<sup>95</sup>**  
Reg. \$62.80

555-8  
4 Litres

**C/D #2 BASEPLATE WAX**



**\$52<sup>50</sup>**

0252  
Pink (5lb)

**C/D FLASKS - UPPER & LOWER**



**\$105**

44AUPPER

**4 + 1 FREE!**

44ALOWER

**C/D MIXING TIPS**



**\$26<sup>95</sup>**

YM100YB  
Pkg(50)

**C/D ARTICULATOR**



**6 + 1  
FREE!**

105067 - High Arch Chrome **\$8<sup>95</sup>**

**C/D RELINE JIG**



**Receive FREE  
Articulator!**

113145 - Chrome-plated **\$42<sup>00</sup>**

**C/D FLASK COMPRESS**



Special  
**\$89<sup>50</sup>**  
Reg. \$99.50

#15  
Compress complete

**C/D IDP INSULATING & DUPLICATING PASTE**

- Helps in the "Break-out" of acrylic denture work
- Coats, protects and safeguards during the divesting process

**\$42<sup>00</sup>**

46525 - Kit - IDP Investing Compound (10oz)  
& Activator (1/2oz)

## PLASTIPAK SHEETS

- Good for trial packing
- Stretches when pulled
- Leaves no wrinkles

**\$15<sup>40</sup>**

45515  
5"x5" Plastic Sheets  
(1000)



## BUFFING AGENT

Designed to give acrylic dentures the highest possible luster with the least amount of time and effort. It is also an excellent polish for gold.



**\$28<sup>40</sup>**

41940

**Buy 2 Receive  
FREE Plastipak!**

## SOLDERING MACHINE

An invaluable time and labor saver plus exceptional laboratory profit generator. Standard equipment for any dental lab.



Special  
**\$985<sup>50</sup>**

Reg. \$1095

42840  
Soldering Machine

## DIAL CALIPER

- Accurate high quality caliper
- Highly accurate to 0.05 mm
- Made of best quality stainless steel



**\$175<sup>00</sup>**

46235  
Dial Caliper

## RAPID - GLAZE DIAMOND POLISH

Eliminates the need for re-glazing porcelain appliances after minor adjustments.



Special  
**\$29<sup>95</sup>**

Reg. \$35.75

69302  
Rapid Glaze (2g)

## WAX IRON

Cast of chrome cobalt alloy. Indestructible in a constant flame.



**\$24<sup>95</sup>** /ea.

**Buy 2 Receive  
FREE Plastipak!**

49102  
Wax Iron Large

Buy 6  
Get 1  
FREE!

# DIATIT TUNGSTEN CARBIDE BURS

Smooth  
Running  
& Longer  
Lasting

Diatit tools feature the advanced Diatit wear protection coating. This is a special material which is added into the surface of the bur after it has been produced. It hardens the tool surface and reduces the surface friction. Results in a smoother running and longer lasting bur.



Material	Tungsten carbide	Diatit
REF	D200 KG 23	<b>\$28<sup>95</sup></b>
ISO-No.	509 104 200220 023	



D200 KG 23 for finishing of CoCr frameworks.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr denture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	D237 KG 65	<b>\$48<sup>50</sup></b>
ISO-No.	509 104 237220 065	



D237 KG 65 excellent removal of material and smooth running for efficient processing of resin.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr denture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	D263 KG 60	<b>\$50<sup>95</sup></b>
ISO-No.	509 104 263220 060	



Due to its shape and the coarse cut the D263 KG 60 is well suited for grinding of plaster dies..

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr denture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	D274 KG 60	<b>\$53<sup>25</sup></b>
ISO-No.	509 104 274200 060	



The D274 KG 60 is suitable for all materials due to this coarse cut.

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr denture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	D292 KG 23	<b>\$28<sup>95</sup></b>
ISO-No.	509 104 292220 023	

Application field	Plaster	Denture resin	Veneer resin	Precious metal/pd-based	CoCr denture/NPM	Ceramic
Working speed rpm	10-20,000	8-12,000	10-17,000	15-20,000	15-20,000	15-20,000



Material	Tungsten carbide	Diatit
REF	D194 KG 23	<b>\$27<sup>60</sup></b>
ISO-No.	509 104 194220 023	

REF	D194 KG 40	<b>\$39<sup>95</sup></b>
ISO-No.	509 104 194220 040	

REF	D194 KG 50	<b>\$49<sup>95</sup></b>
ISO-No.	509 104 194220 050	



D194 KG 23 for efficient CoCr processing



The coarse cuts of the D194 KG 40 (picture 3) and the D194 KG 50 (picture 2) guarantee fast and perfect finishing of resin.

# Multisil-Soft

The permanently soft relining system.

**NEW!**

bredent group

REF 540 0104 5

**SAVE \$40**



- Assortment  
 Multisil-Soft  
 2 x 50 ml Multisil-Soft in cartridges  
 5 ml Multisil-Primer  
 10 ml Multisil sealing liquid  
 12 pieces Mixing cannulas  
 1 piece Silicone burs

Special Intro Price  
**\$225**

Regular Price \$265.00



reliable bonding



ready-to-use



grindable

# DTK Adhesive Kit

This self curing adhesive is ideal for joining attachments and is excellent for resin bonding metal to metal. The paste / paste formulation provides a consistent mix and a high bond strength. The anaerobic setting properties also provide an extended working time. Licensed under Health Canada Regulations to provide you with the peace of mind that you and patient deserve!

Regular Price \$166.00  
 Special  
**\$126**

**SAVE \$40**



Adhesive Kit - 54000106

# Smart Pins

Narrow, short model pin with optimum insertion into metal sleeve.

- precise, conical fit of the pins into the sleeve.
- produces flat bases.
- reduced plaster processing ensures reduced expansion within the model.
- ideal for crowded spaces e.g. anterior area.
- particularly narrow retention head [1.6 mm] for delicate dies.



Special  
**\$75** Reg. \$91.50

Pins: 366-2000  
 Pkg. [1000]


Special  
**\$79<sup>95</sup>** Reg. \$95.55

Sleeves: 366-2100  
 Pkg. [1000]


# NEW! Renfert Polish

## Diamond ceramic polishing pastes

Renfert Polish all-in-one: for polishing veneering porcelain, hybrid ceramics, high performance resin; as well as for zirconium oxide and lithium disilicate ■ Renfert Polish  $ZrO_2$ : diamond granulation especially suited for zirconium oxide ■ Renfert Polish  $LiSi_2$ : diamond granulation especially suited for lithium-disilicate



Renfert Polish  $ZrO_2$ ,  
10 g  
Item Code 510 1000  
Special price:  
**\$97<sup>25</sup>**



Renfert Polish  $LiSi_2$ ,  
10 g  
Item Code 510 2000  
Special price:  
**\$97<sup>25</sup>**



Renfert Polish  
all-in-one, 10 g  
Item Code 510 0000  
Special price:  
**\$97<sup>25</sup>**

**Just Prestin  
CDT and laboratory  
owner:**

“With the Renfert Polish all-in-one you reach the desired gloss level in a few seconds – the best all-round polishing paste I know!”

# NEW!

## GEO Expert Wax Set

### Design Wax Set nach der Philosophie von August Bruguera

Easy reproduction of natural dental forms and colors thanks to the how-to-guide by August Bruguera ■ Quick and clean fabrication of dental wax-ups due to excellent modeling and scraping characteristics of modeling and effect waxes ■ Precise building up even of finest structures such as cusps and mamelons thanks to excellent flow and solidification properties and optimum stability



GEO Expert Wax Set,  
Item Code 634 0001  
Special price:

**\$155<sup>15</sup>**

*A. Bruguera*



**August Bruguera  
Opinion Leader**

»With the GEO Expert Wax Set that I co-developed, you receive an exhaustive how-to-guide intended for experts, individualists and for all those who are willing to learn. I show the monolithic, diagnostic, hyperreal/natural waxing-up, step-by-step and in detail, clearly illustrated by photographs and texts.«



# Get Free MLB jersey of your favourite team and players!



Purchase a..

- Silent TS2, Get 3 Free Jerseys!**
- Silent TS, Get 2 Free Jerseys!**
- Silent Compact, Get 1 Free Jersey!**
- Silent CompactCAM, Get 1 Free Jersey!**
- Twister Mixer, Get 2 Free Jerseys!**
- Bacic ECO, Get 1 Free Jersey!**
- MT3, Get 1 Free Jersey!**
- Top Spin, Get 1 Free Jersey!**



..all team Jerseys can be ordered, with name & number of your favorite player!

## Twister / Twister evolution

Vacuum mixing unit

The high torque of the mixing motor ensures reliable processing of large amounts of mixing material. Intuitive programming and operation using a large, highly luminous display.



**\$2730** Twister Evolution  
Item Code 1828 1000

**\$2150** Twister Mixer  
Item Code 1826 1000

## **NEW!** SILENT compact

Dust extractor for single workstation

No need to change the dust bag as fine filter technology takes place without a filter bag and thanks to mechanic filter cleaning. Appliance controlled via an intelligent, automatic switch. Noise level reduced due to a fully encapsulated turbine unit.



SILENT compact,  
Item Code 2934 1000  
Special price:

**\$1458**

## SILENT TS2

Two workstation extraction

Suction power only where and when it is needed, thanks to the intelligent independent control of the 2 separate pinch valves. Constant suction power thanks to automatic re-adjustment of motor output. Intuitive, convenient operation with large, clear multifunctional display.



SILENT TS2,  
Item Code 2930 1050  
Special price:

**\$3985**

## Top spin

Pin drilling unit

Unique shaft bearing technology ensures a drilling accuracy of  $\leq 0.01$  mm. Powerful motor without losing performance during drilling 8,000 rpm. Unique drill technology enables extremely easy drilling without resistance.



Top spin,  
Item Code 1840 1000  
Special price:

**\$1528**

## **NEW!** Basic eco

Fine sandblasting unit with 2 tanks

Maximum field of vision thanks to a large viewing panel with innovative LED lighting. Process reliability based on many years of experience in dental sandblasting technology. Saving in abrasives using optimal air-abrasive mixture.



Basic eco,  
Item Code 2949 4025  
1 x 25-70 µm, 1 x 70-250 µm  
Special price:

**\$980**

## MT3 Klettfix

Wet trimmer

Effortless trimming due to maximum motor performance. Safety: when the front door is opened a safety switch stops the motor and water. Uniform coverage by the optimally positioned spray tube prevents clogging of the trimmer disc.



MT3 incl. Klettfix  
trimmer disc,  
Item Code 1808 1000  
Special price:

**\$1275**

# 8000CD - AUTOMATIC STEAM CLEANER



The newest member to our professional series steam boilers, the 8000CD provides continuous steam for an unlimited amount of time as it must be hooked up to a direct water line. There's no more down time waiting for the boiler to re-pressurize. The 8000CD is the ideal choice for environments that require large volumes of uninterrupted steam.

## WHY CLEAN WITH STEAM?

The only ingredient required for the 8000CD is water. Cleaning has never been safer, healthier or more environmentally friendly.



**\$1799**

8000CD



# AcryLux

*We take care of your smile*

- Best value in the market for Three-Layer Tooth.
- Large selection of moulds (37 Anterior & 8 Posterior & Shades including Bleach).
- New everyday Price of only **\$10.50** per Card.

3D View dentin/neck



3D View enamel/neck




3D section of stratification



**SPECIAL OFFER!**

Live Working Mould Guide 37 Anterior & 16 Posterior only

**\$499**



3 LAYERS OF ACRYLIC





## HAN-D-III WAXER



01420A

Special

**\$175<sup>00</sup>**

Reg. \$190.00

## STAIN MIXING PLATE



**\$27<sup>95</sup>**

04189-19

## LABO-MATE 80



**6 + 1  
FREE!**

**\$42<sup>95</sup>**

02030



## SONIC DENTURE CLEANER




**\$8<sup>95</sup>**

04210

**12 + 1  
FREE!**

## POLISHING PATS



04119  
Super Hard



04120  
Hard



04121  
Regular



04118  
Fine

**\$27<sup>50</sup>**

## MICRO-TORCH



**\$31<sup>95</sup>**

04020



## NITRILE POWDER-FREE EXAM GLOVES

**\$16<sup>95</sup>** / Box



CODE/SIZE
92395-XS
92396-S
92397-M
92398-L
92399-XL

**200**  
Gloves / Box

*Double the Normal Amount! (100/box).*

**9 + 1 FREE!**

## CHROMA 2000 RETAINER BOXES

Simple, economical, durable retainer boxes. Hinge design with snap lid.

**4 + 1 FREE!**



CR20002X Blue (12/box)  
CR20005X Salmon (12/box)



“The **BEST** there is!”  
Laminated Mouthguards



## PROFORM MACHINE IV Vacuum Former

#7000338

**\$599<sup>00</sup>**

## PROFORM Mouthguard Material

Includes: Solid, 2-Colour & 3-Colour (12/box)

**4 + 1 FREE!**



- ✓ Anterior brace for extra protection!
- ✓ Custom-fitted for great oxygen intake!
- ✓ Easy to attach break-away strap!
- ✓ Custom-made shock-absorbing surface!
- ✓ Professionally made for a perfect fit.
- ✓ Laminated material for superior strength.
- ✓ Excellent abrasion resistance!
- ✓ Made by dentists worldwide.

# ZHERMACK Elite HD



Buy 8 Packages  
and Receive 1  
Package (50)  
Mixing Tips **FREE!**

**4 + 1  
FREE!**



# ZHERMACK Alginates

**5 + 1  
FREE!**

Purchase 20 and  
Receive a **FREE**  
Bowl and Spatula!



# ZHERMACK Elite Double

**2 + 1  
FREE!**

Buy 2 and Receive a **FREE**  
Duplicating Flask!



020186



**\$30  
Value!**

**ANY SIZE!**

# core3d centres

## Core3d Milling Centres

– It's about connections.

This past summer Central Dental Ltd. and Core3d Milling Centres entered into an exclusive sales partnership for Canada so that we could provide access for dental technicians and denturists to the means to respond to the avalanche of change in dental technology being driven by digital technologies.

Every day more dentists are demanding prosthetic solutions that are digitally driven, from full-contour zirconia crowns to milled implant abutments to printed models for C&B work derived from intra-oral scanners. Accessing the workflows and processes to obtain these solutions for your clients can often be confusing and intimidating. Do you have the right scanner and software? Can the digital files you receive or create be received and processed by the outsource partners you would like to work with? What product options do you have to select from? Are you getting the level of technical support needed to successfully complete your more complex cases?

Core3d Milling Centres have the processes and systems that allow you easily and efficiently connect to a wide range of solutions, and are committed to levels of technical support and product quality that will allow to not only meet the demands of your existing clients but also reach out into your market place to connect with new dentists/patients demanding digitally-driven dentistry.

### Some of the Core3d solutions are:

- Full contour zirconia crowns and bridges
- Emax crowns and bridges
- Custom-milled implant abutments
- Custom-milled implant bars
- Milled models for Itero scanners
- Printed models for a range of intra-oral scanners

Ask us about  
cubeX<sup>2</sup> Super  
High Translucent  
Zirconia!

**DDcubeX<sup>2</sup>**

For help in understanding and connecting to the full range of solutions from Core3d please speak to your Central Dental Ltd. sales representative or call us at -1-800-268-4442.