



CLINICAL SYSTEMS

STR SERIES

STANDARD FEATURES

- 0.90mm CNC manufactured DKP steel
- Electrostatic powder coating
- Various RAL color range
- 30mm laminate post-formed worktop with 3mm PVC side cover
- Modular system, compact design
- Soft close drawers, plastic trays
- Waste bins on sliding tray
- Autoclave sections
- Upper cabinets and dispenser cabinet
- Long lever faucet, acrylic sink
- Chrome accessories
- CE certified

OPTIONAL FEATURES

- Compact, wood (beech), acrylic and stainless steel worktops
- Ceramic, inox and glass sinks
- Infra-red water tap
- Wall panel



STR - 5

60 x 123 x 205cm (h)
Weight: 140 kg

\$2495.00



\$1195.00

STR - COLUMN

60 x 60 x 205cm (h)
Weight: 55 kg

AQUA SERIES

STANDARD FEATURES

- 0.90mm CNC manufactured DKP steel
- Electrostatic powder coating
- Various RAL colour range
- Modular system, compact design
- 19mm thickness, paintable glass, built-in sink
- Soft close drawers
- Waste bin inside sink cabinet
- Long lever chrome faucet and chrome pop up siphon
- Complete chrome accessories
- Plastic Trays
- CE certified



AQUA - 3

56 X 166 X 89cm (h)
19mm glass, built-in sink
Waste bin
Long lever faucet & pop-up siphon
Weight: 150 kg

\$3895.00



AQUA - 2

56 X 117 X 89cm (h)
19mm glass, built-in sink
Waste bin
Long lever faucet & pop-up siphon
Weight: 100 kg

\$2795.00

AD7 Mini Sterilizer

The most compact and convenient autoclave for surgical instruments, high-speed handpiece and trays

- One-click-ease to enable 4 pre-programmed sterilization cycles from the user-friendly control panel
- Speedy 3.5 minutes sterilization cycle available for urgent usage
- Automatic door opening mechanism
- Hassle-free auto water drainage design
- 6 safety mechanisms for utmost operation safety
- High quality stainless steel chamber for longer reliability
- 14 handpieces can be disinfected effectively in every cycle.
- Two containers (HP rack, HP basket) are provided as standard accessories
- Dimensions: (L) 37.5 cm x (W) 32 cm x (H) 34.5 cm
- Chamber: (D) 12.5 cm x (L) 24 cm



\$2495.00

code: AD7MINI

The *Miracle* of New Ultrathin Nitrile!

NITRILE POWDER-FREE EXAM GLOVES

A truly superior performance you would never expect from a nitrile glove! Optimal stretchability, just like latex! Long-lasting comfort, durability and decreased hand fatigue. Beaded cuff snugs wrist, eliminates roll-up, minimizes entry of liquids. Thickness: 2.5 mil. Finger-textured finish for precision gripping. Latex-free, non-allergenic.

Color: Teal Green.

CODE/SIZE
92395-XS
92396-S
92397-M
92398-L
92399-XL



Intro Special

\$16.95 / Box



SR Ivocap Hi Impact

Acrylic Capsules
Shades Available:

Preference - 559615AN

US-L - 531402AN

US-P - 531403AN

(Pkg. 50)

Special

\$349.95

Reg. \$415.00

SAVE
\$65





Gradia Gum



Buy Starter Kit & Receive
1 Primer FREE!

\$80 VALUE

Pattern Resin



Buy 1 Pattern
Resin Powder
Receive 1 Liquid
FREE!

Initial Ti

Buy Initial Ti Basic Set
Receive Initial Ti
Bonder FREE!



FujiRock Optiscan



**NEW Scannable
Die Stone!**

Buy 2 FujiRock Optiscan
5kg Receive 1 FREE!

Nature-Cryl Hi Impact

Buy 25 lbs Natural Super Hi
Impact P&L Receive Gradia
Gum Starter Kit FREE!
Or Receive FREE Liquid!



OR

Buy 5 lbs. Natural Super Hi
Impact Powder Receive
1 qt Liquid FREE!

Nature-Cryl MC



Buy 5 Lbs of Nature-Cryl MC
& Get 1 qt Liquid FREE!

Fuji Rock EP

12 Kg.



3+1

890226 - Brown
890227 - Pastel Yellow

Please fax your invoices to GC to receive free product: 1 (800) 423-2963



Zetalabor & Titanium Putty



BUY 3 + 1 FREE!

ANY SIZE!

Hi Tray's BUY 1 + 1 FREE!



Elite HD

BUY 4 + 1 FREE!



Alginates

BUY 5 + 1 FREE!



Elite Double BUY 2 + 1 FREE!

ANY SIZE!



Occlufast BUY 2 + 1 FREE!



Gingifast BUY 2 + 1 FREE!



Platinum Putty BUY 2 + 1 FREE!



HOLIDAY SPECIALS!



• **FiBER FORCE CST**®



Inspired by SIMPLICITY AND AFFORDABILITY

The CST® technique makes the fabrication process both fast and simple. A CST® framework can be made in approximately 30 minutes and requires no casting, scanning, or milling – making it accessible to any dental professional with a light curing unit. All the usual protocols for fixed-hybrid dentures are followed which means familiar routines and processes are maintained.

Not only are CST® frameworks strong and easy to fabricate, they're affordable – making a fixed denture solution more accessible to more people.

FREE pack of mesh, UD, 2mm or 1mm Fiber Force with the purchase of a CST intro kit. Promo will be in effect until Dec 31st.

CST® complete kit (CST50): \$475

Contents: 2 - 17.75" (450mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:6), 2 - 11.75" (300mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:4), 1 - 3ml syringe of FIBER FORCE® CST® Flow Pink Resin with tips, 1 - 5ml bottle of FIBER FORCE® CST® Bond Resin, 4 - FIBER FORCE® CST® Temporary Support Pillars & Instructions.

CST® refill kit (CST60): \$175

Contents: 1 - 17.75" (450mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:6), 1 - 11.75" (300mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:4), 2 - FIBER FORCE® CST® Temporary Support Pillars & Instructions.



To Order Call: 1.800.250.5111

Qu system



Bredent's Qu-System

Qu-Resin: What everyone is talking about!

- Pick up attachments intraorally.
- Verification of splints.
- Utilize for open tray pick-up on implants.
- Excellent denture repair resin that bonds perfectly between the denture base and high impact acrylics using the Qu-Resin connector.
- Two component system is based on a diacrylate and sets in 3 minutes.
- Repairs can be done chairside.
- Available in pink or dentin to accommodate repairs of the denture base or teeth!

- (54001161) Qu-resin pink kit **\$125.95**
- (54001162) Qu-resin dentin kit **\$125.95**
- (54001165) Qu-resin pink, 50 ml **\$99.00**
- (54001166) Qu-resin dentin, 50 ml **\$99.00**

Bredent's Pi-Ku-Plast hp 36

Perfect for highly precise implant work, Pi-Ku-Plast 36 modelling resin has the lowest contraction of any pattern resin in the market with a shrinkage value of only .036%. It also provides optimum contouring properties since it will not slump. A clean burnout provides top quality casting results.



- (54000219) Blue
- (54000220) Red
- (54000216) Transparent *
- (54000217) Yellow *
- (54000218) Orange*

Special
\$259.95

Reg. \$287.61



** Yellow, Orange and Transparent are special orders.*

Bredent Exacto-Form: Polyethylene modelling resin

The perfect base for precision crown and bridge construction. Exakto-Form's minimal shrinkage ensures a high accuracy of fit. The excellent flowability avoids bubbles and forms accurate margins.

Available in five colours:

Isolating liquid CODE: (52000210)
\$31.30 / 125 ml



- green (52000280)
- blue (52000281)
- ivory (52000282)
- yellow (52000284)
- grey (52000283)

Special
\$99.95

kit contains 6 tins each
of component A & B



Abraso-Sil Acryl

Ø 80 mm
1 piece

This buff consists of a nonwoven fibre fabric in the centre between two layers of silicone-coated cotton fabric. On the outside there are two rows of bleached Chungkink bristles. This brush absorbs polishing paste or pumice and releases it slowly for efficient prepolishing.



Special
\$14.95
Reg. \$21.22

code: 350 0099 3

Set-up grinding tool

1 piece

SAVE \$30

- Two grinding tools in one.
Grinding without exchanging tools in a single working step.
- Quick adaptation of the underside of the tooth to be set up.
 - Grinding in of occlusal stops.

code: 340 0101 0

Special
\$59.95
Reg. \$89.39

occlusal



The small, precisely shaped grinding tip with fine, perfectly cutting diamond grains provides the ideal precondition for well-aimed and rapid grinding in of occlusal contacts.

basal



The large grinding area with its optimized shape and selected natural abrasive diamonds ensures maximum removal of material and thus accurate and quick grinding.

Bredent Diagen Turbos - Fine Grit!

Diagen turbos are very aggressive sintered diamonds that work effectively on zirconium oxides such as Procera™, Inceram™ and all other ceramic materials. Each stone is moulded to provide optimum grinding properties and cuts effectively down to the last layer.



34000150	34000155	34000160	34000165	34000170	34000180	34000190	34000210	34000220	34000240	34000250	
\$60.00/2pk	\$60.00/2pk	\$68.00/2pk	\$79.00/2pk	\$85.00/2pk	\$85.00/2pk	\$85.00/ea.	\$64.00/ea.	\$64.00/ea.	\$65.00/ea.	\$63.00/ea.	



HOLIDAY SPECIALS!



DENTSPLY'S Q4 BUNDLE #1

PURCHASE \$9,000 IN DENTSPLY REMOVABLE PRODUCTS, GET A MACBOOK AIR 11" 128GB (VALUED AT \$999) FREE!

PDR13-Q4-9



DENTSPLY'S Q4 BUNDLE #2

PURCHASE \$4,500 IN DENTSPLY REMOVABLE PRODUCTS, GET AN iPad4 WITH RETINA DISPLAY 16GB (VALUED AT \$500) OR A MICROSOFT SURFACE TABLET RT WITH KEY BOARD (VALUED AT \$450) FREE!

PDR13-Q4-10

DENTSPLY BUNDLES

Eligible products: GENIOS & Portrait IPN Denture Teeth, Lucitone 199, Triad, Al-Cote, Luci-Sof, Lucitone FRS, Eclipse, Neoloy, Enterra, Visco-Gel, Success, INTerra, Faskut, TurboMAX, Fas-Por+, Hy-Pro, PermaSoft. Visco Gel and Druformat Scan, Portrait IPN and GENIOS Living Mould Guides. For redemption, please fax invoice to DENTSPLY at 905-851-9809.



DENTSPLY'S Q4 BUNDLE #3

PURCHASE \$2,000 IN DENTSPLY REMOVABLE PRODUCTS, GET A 4TH GENERATION iPod TOUCH 16GB (VALUED AT \$229) FREE!

PDR13-Q4-11



DENTSPLY'S Q4 BUNDLE #4

PURCHASE \$1,100 IN DENTSPLY REMOVABLE PRODUCTS, GET A CHOICE OF A FREE:

- LEAPPAD2 EXPLORER LEARNING TABLET OR
- 5X TIMOTHY'S K-CUPS REFILLS OR
(MILD ROAST BREAKFAST BLEND, 24-CUPS/PKG)
- KEURIG ELITE COFFEE MACHINE

PDR13-Q4-12

Dealer Promotions: All Free Goods which are gift items will be shipped by DENTSPLY Canada to the End User upon receipt of invoice copy as proof of purchase. Invoice must be dated October 1st, 2013 to December 31st, 2013. Invoice copies are to be mailed directly to DENTSPLY Canada Attn: Marketing Department or faxed to 905-851-9809 or email promotions.canada@dentsply.com. **Invoice copies must be received no later than January 31st, 2014 to qualify.** Purchases made as related to free gift items are strictly not returnable for credit. Unless indicated as shipped up front to dealer, all free goods which are product related are to be shipped by the dealer and claimed back from DENTSPLY Canada with proof of sale to the End User. Purchases for the above promotions do not qualify for any other DENTSPLY special programs. Offers are valid in Canada only.

DENTSPLY.CA

HOLIDAY SPECIALS!



GENIOS MOULD GUIDE (ANTERIOURS - #01301ZAANT)

Regular retail price: \$987.52
Promotional price: \$550.00
SAVINGS OF \$437.52 (44%)

Contains: A Genios shade Guide with 16 shade tabs, 15 anterior moulds (9 upper & 6 lower) in A1 and B1 shades. Plus: Two mould Charts (English and French). By ordering through this special promotion, you can customize your GENIOS Mould guide with the Shades and Moulds of your preference.

PDL13-Q4-2



GENIOS MOULD GUIDE (ANTERIOURS & POSTERIOURS - #01301SM) CUSTOMIZABLE

OFFER FOR NEW USERS!
Regular retail price: \$1,483.74
Promotional price: \$999.00
SAVINGS OF \$484.74 (32%)

Contains: A Genios shade Guide with 16 shade tabs, 18 posterior moulds (9 upper & 9 lower) plus 15 anterior moulds (9 upper & 6 lower) in A1 and B1 shades. Plus: Two mould Charts (English and French). By ordering through this special promotion, you can customize your GENIOS Mould guide with the Shades and Moulds of your preference.

PDL13-Q4-3



PORTRAIT IPN MOULD GUIDE (#902745C)

OFFER FOR NEW USERS!

Regular retail price: \$1700.00
Promotional price: \$999.99
SAVINGS OF \$700.00 (37%)

Contains: Anteriors Upper - 11H, 12E, 12G, 21D, 21X, 22E, 32E, 42F, 45H, Anteriors Lower - C, E, F, H, K1, N, P, R, S, Posteriors 20 Deg (upper and lower) - 29L, 29M, 31L, 31M, 33L, 33M, 35M, Posteriors Biostabil 22 Deg (upper and lower) - 530, 532, 533, Posteriors Anatoline (upper and lower) - 632, Posteriors 0 Deg (upper and lower) - 632, Portrait IPN Shade Guide with 16 A-D, 3 bleached & the top 8 Bioform shades. Available in A1, A2 and B1 shades.

PDL13-Q4-4



DENTSPLY REPAIR MATERIAL

Purchase one DENTSPLY Repair Material Labpak Kit*, **RECEIVE ONE 2-LB (900G) TUB OF TRIXA PUTTY AND ONE 41G TUBE OF TRIXA CATALYST (#905449) FREE! A SAVINGS OF \$74.55!**

*Shades: Pink #680010, Pink Fibered #680011, Light Fibered #680028. For a demo on Repair Material, visit DENTSPLY Prosthetics' YouTube Channel

PDL13-Q4-5



LUCITONE FRS

Purchase two 12-packs of Lucitone FRS*, **RECEIVE ONE BEVELED FASKUT LARGE ABRASIVE WHEEL (#904557) AND ONE FASKUT ALUMINUM OXIDE ABRASIVE POINT ASSORTMENT (#99116A) FREE! A SAVINGS OF \$87.89!**

*Shades: Original #904854, Light Pink #904855, Light Reddish Pink #904856, Dark Pink #904857, Clear #904853

PDL13-Q4-6



TURBOMAX ALGINATE AUTO MIXER

Purchase a TurboMAX Alginate Auto Mixer, **RECEIVE 2 CANS OF JELTRATE CHROMA OR JELTRATE PLUS FREE!**



DENTSPLY.CA

TECHNICAL CORNER

Chairside & Laboratory repair techniques using Bredent Qu-Resin

Qu-Resin is a Quick Setting, Self-Curing denture repair resin based on diacrylate for Intraoral and Extraoral (Direct and Indirect) use.

LABORATORY TECHNIQUE



1
Tooth no. 24 is extracted and the partial denture is "transformed" into a full denture.



2
A key is prepared to record the situation.



3
The areas to be completed are roughened, Qu-connector is applied and light-cured for 90 seconds. Then Qu-resin dentin is injected into the key; the key is replaced again and the tooth shape is completed from the oral direction.



4
After curing, remove the key and complete the gingival sections of the denture with Qu-resin.

- Available in 2 colours (dentin / pink)
- Sets in 3 minutes
- For Intra-and extraoral (direct and indirect use)
- Self-curing resin in the cartridge, no pressure pot required
- Connector for bonding with other resins; can be activated with all UV units
- Excellent processing and polishing properties

CHAIRSIDE TECHNIQUE

Qu-resin for intraoral use.
 The dentist can perform chairside repair of the denture and reduce the amount of time considerably.




Lower CoCr restoration is "transformed" into a temporary full denture after extracting the abutment teeth.



Impression of the actual situation with alginate.



After the extraction of the teeth, the positions of the missing teeth in the impression are filled with Qu-resin.



The areas to be completed in the transition zones are conditioned with Qu-Connector.



Missing pink-colored areas are filled directly in situ.



Material can be directly injected below small, non-supported areas in the mouth.



Denture completed with Qu-resin dentin and Qu-resin pink.



Completed temporary denture in situ.

Plastic Cups



#30-5148

5 oz Cup
(Package of 1000)

Special
\$28.75
Reg. \$31.95



Headrest Covers

#24-1114 (Case of 500)



Special
\$28.80
Reg. \$32.00

#24-1114 (Headrest Cover Tissue/Poly -10"x13")

Dental Bibs

- #23-1210 - Blue
- #23-1211 - Green
- #23-1212 - White
- #23-1213 - Rose
- #23-1214 - Yellow



Special
\$19.95 (Package of 500)
Reg. \$23.95

Sterilization Pouches

Special
\$14.65
Reg. \$16.30
#P0140255



Special
\$8.55
Reg. \$9.50
#P090135

#P140255 - 5.25" x 10" Package of 200
#P090135 - 3.25" x 5.25" Package of 200

Nov 1-2, 2013
DTABC Annual Convention
Sheraton Vancouver Airport Hotel
Richmond, BC

Nov 7, 2013
FiberForce CST
Cable Stayed Technology
Hands-on Course
Scarborough, ON

Nov 9, 2013
DTA Lecture
Dental Technology
Academy Lecture: Business
Scarborough, ON
@ Central Dental Ltd.

Nov 15, 2013
Phonares II Setup Course
(9:00am - 5:00pm)
Denture Setup using Phonares II
denture teeth - Hands-on Course
Scarborough, ON
@ Central Dental Ltd.

Nov 16, 2013
DISC Lecture
Denturism Implant Study Club
www.disccanada.ca
Scarborough, ON
@ Central Dental Ltd.

Nov 28, 2013
FiberForce CST
Cable Stayed Technology
Hands-on Course
Vancouver, BC

Nov 30, 2013
FiberForce
Fiber Framework System
Hands-on Course
Vancouver, BC

Holiday Hours of Operation
Closed on: Dec 24, 25, 26, 31 & Jan 1

For upcoming events and more please visit our
website at www.CentralDentalLtd.com