




The visio.lign veneering system is based on multi-layer veneers (novo.lign) and developed in accordance with the example of natural strength. A bonding system with matched shades and suitable individualization, and add-on materials enables the perfect combination of individual and efficient veneering techniques. With abrasion values comparable of natural teeth, visio.lign provides gentle occlusion - an advantage especially for implant-supported restorations.



Visio.lign Hands on Course Production of a 4-unit Anterior Bridge on Zirconia

- Introduction to Visio.lign
- Description of the Visio.lign Components
- Case Planning (Short Break)
- Set-up frameworks using high Impact PMMA Visio.lign veneer model
- Framework Design and Conditioning using COMBO.LIGN Opaquer
- Conditioning of veneers using VISIO.LINK
- Luting veneers with COMBO.LIGN dual-hardening adhesive composite (Lunch Break).
- Final shaping and dentitive design using the CREA.LIGN nanolled composite including pink custom gingiva (Red & White Aesthetic)
- Polish frameworks using the VISIO.LIGN TOOL KIT

Friday, October 10th
Toronto ON

Friday, November 7th
Toronto ON

Thursday, November 20th
Calgary Alberta

Friday, November 21st
Vancouver BC

TALLADIUM

MASTERmill Dental Milling Machine with Dental Wings Scanner



Key Features & Benefits:

- Zirconia, pmma and wax dental milling machine
- Mill comes with laptop and cad/cam software included
- Easy to operate with user-friendly dentalcam software
- No previous milling experience required
- Extremely compact desktop design
- Automatic tool changer with tool length detector
- In-built monitoring systems: tool breakage monitoring, compressed air monitoring.
- The jäger high frequency spindle (i.E. The most important part of the machine) is isolated from the cutting chamber. (This prevents the abrasive zirconia dust from damaging the spindle)
- Quiet, precise, & clean

Price
\$19,000

Price
\$31,500

Lease from:
\$1008/month
5 year term

Talladium Luminesse HT 95mm

5170-HT	95mm X 10mm	\$155
5171-HT	95mm X 12mm	\$178
5164-HT	95mm X 14mm	\$210
5165-HT	95mm X 16mm	\$240
5172-HT	95mm X 22mm	\$285

Talladium Luminesse HT 98mm

5142-HT	98mm X 10mm	\$164
5143	98mm X 12mm	\$182
5144	98mm X 14mm	\$214
5145	98mm X 16mm	\$243
5146	98mm X 18mm	\$250
5147-HT	98mm X 20mm	\$267
5148-HT	98mm X 25mm	\$318

Talladium Luminesse HT 100mm

5139	100mm X 10mm	\$164
5135-HT	100mm X 12mm	\$190
5136	100mm X 14mm	\$224
5137-HT	100mm X 16mm	\$251
5140	100mm X 18mm	\$259
5141	100mm X 25mm	\$277

Luminesse HT CAD/CAM Discs



Talladium is moving forward in the milling era!

Our yttria-stabilized zirconia is manufactured using the highest quality powder offering amazing translucency like no other! Our blocks and discs are isostatically cold pressed and partially sintered producing the highest quality fit available for single units up to full frameworks with absolutely no distortion. Talladium offers 18 different blocks and 13 different discs to choose from. Available for the Sirono Inlab, KaVo Everest, Weiland, DentMill and Zirconc Zahn machine.



5198	98mm x 10mm	\$57.50
5200	98mm x 12mm	\$57.50
5201	98mm x 14mm	\$57.50
5202	98mm x 16mm	\$57.50
5203	98mm x 18mm	\$57.50
5204	98mm x 20mm	\$57.50
5205	98mm x 25mm	\$76.70

5199-100	100mm x 10mm	\$57.50
5206-100	100mm x 12mm	\$57.50
5207-100	100mm x 14mm	\$57.50
5208-100	100mm x 16mm	\$57.50
5209-100	100mm x 18mm	\$57.50
5210-100	100mm x 20mm	\$57.50
5211-100	100mm x 25mm	\$76.70

Luminesse CAD/CAM Wax Discs

- Talladium, Inc.'s Luminesse Green CAD/CAM Wax Discs are strong and stable.
- Offer a harder consistency to create high accuracy and precision fit margins for copings and full crowns.

3/pk

5175 - Diagnostic White Wax
100mm x 20mm for milling (1/pk)
\$34.50



Waxlectric II

Waxing-up is simple, precise and stress-free due to the accurately set, constant and evenly distributed temperature of the instrument tip. Huge saving in time because separate heating of the instrument and material is no longer required.

Special
\$399⁰⁰

Reg. \$476.46

2157-1000

Waxprofi

The Waxprofi produces improved prosthetic results more efficiently. Instruments and wax are no longer heated separately, resulting in distortion-free wax denture bases and a huge saving in time.

Special
\$329⁰⁰

Reg. \$347.30

1440-1000

Twister

The compact Twister vacuum mixer can be employed for all conventional casting and modelling materials used in dental technology. All mixing parameters are quickly and easily input. The high performance motor and high quality vacuum pump guarantee homogenous and bubble-free mixing results for dental plasters, investments, and silicones.

Special
\$1775⁰⁰

Reg. \$1879.40

1826-0000

Basic Classic

Fine sandblasting unit with up to 2 tanks.

Special
\$1025⁰⁰

Reg. \$1141.68

2947-4025

MT3

Wet trimmer

1808-1000

Special
\$925⁰⁰

Reg. \$975.00

Dustex Master

Dust and vapours are extracted directly where they originate: at the restoration. The tilt and slip resistant Dustex master plus protects the technician from splinters, dust and vapours. The sophisticated, ergonomic design and the armrests provide a comfortable working posture for fatigue-free work.

Special
\$580⁰⁰

Reg. \$613.55

2626-1105

Hotty LED

Wax dipping pots with a controllable working temperature for gently heating dipping waxes. High-grade wax components are heated without impairing the quality.

Special
\$285⁰⁰

Reg. \$306.60

1461-1000

Sympro

This cleaning device is a new development in the well-known pin cleaning technique. Its extraordinary efficiency results from the blend of the unit's characteristics, the specifically developed abrasive materials and the cleaning fluids. The fully automatic and highly efficient cleaning process doesn't take up any of your valuable resources for unpleasant tasks, that can however be offered as a value added service.

Special
\$895⁰⁰

Reg. \$1027.00

6500-0000

BEAT THE FREEZE!

Galaxy II™ Universal Investment

Superior Quality - Flushed with nitrogen for ultimate freshness. Ceramic fillers provide ultra smoothness.

Reliability - Precise expansion for a plug-in fit - every time.

Consistency - Every batch is CDT approved. 28 quality control checks insure a perfect fit.

Versatility - Use small to large rings with confidence. The additional expansion you need for difficult alloys.



Special
\$141.60
Reg. 161.60

**Save \$20
a Box!**

1400-2 - (Package of 144 x 60g. Pouches) or
1400-2X - (Package of 100 x 100g. Pouches)

Diashine

A compound which will give a high gloss on a variety of surfaces, including ceramics, composites and metals. It is manufactured in a range of abrasive grit sizes and is contained in different viscosity carriers, or bases.



FINE-3G

Price
\$59.00



Boss Pillow Box

New crown box made of clear plastic with two holding clear membranes. When a single crown, or a threeunit bridge or a six-unit bridge or more is placed in the box, it appears as floating in the air, even as it is held firmly in place by the membranes. Perfect for shipping all Ceramic Crowns or Porcelain Veneers.



\$1.25

LARGE 3" x 3"

8070-04

.80¢

SMALL 2" x 2"

8070-02

Micro-Kleen3™

Accelerated Kill-Time Formula™



SAVE WITH MICRO-KLEEN³ INTRO KIT

Intro Kit Includes:

- > 1 - 3.78 L pour bottle
- > 2 - 709 mL spray bottles
- > 2 - Canisters of wipes
- > 1 - Coupon for future purchases



Intro Kit Regular
Value: **\$92.00**

Intro Kit Promo
Price: **\$79.95**

Z E R O C O M P R O M I S E

- 1 Most Effective Cleaner**
The most important first step!
If you can't clean it you can't disinfect it!
- 2 True Broad Spectrum Disinfectant**
MK3 is bactericidal, virucidal, fungicidal & tuberculocidal, because you don't know what's on your surfaces
- 3 Short Contact Time without Sacrificing Efficacy**
MK3 will kill 100% of all labeled pathogens in 3 minutes!
- 4 Surface & Equipment Compatibility**
MK3 is safe for use on all healthcare environmental surfaces!
- 5 Safety of Use**
MK3 is safe for use on inanimate surfaces and leaves no hazardous chemical residue when used as directed!

FREE GIFT!

Included with Purchase of the Micro-Kleen³ Intro Kit
Order Code: 038-5997



Receive a
FREE glass nail file with
SWAROVSKI® elements. One gift per
Micro-Kleen³ Intro Kit purchase.

Included in Intro Kit. Retail value: \$15.95 ea.

GC Initial Ceramics

GC Initial can create metal-ceramic restorations all the way to full-ceramic. This product consists of six ceramics in one integrated product line with one color system. No longer does the ceramist have to spend valuable time learning to use a myriad of different ceramics from an array of manufacturers to create the variety of restorations which may be needed to complete a case. Each GC Initial ceramic is adapted to meet the needs of its particular fabrication process and framework, such as its coefficient of thermal expansion (or CTE). Stress cracks and fractures are therefore avoided, which assures strong bonding, adhesion and retention.



4+1 FREE

For Any GC Initial Refill Porcelain

NEW FUJIROCK

One of the world's most popular Type IV Dental Die Stone manufactures now brings you NEW GC FUJIROCK, NEW GC FUJIROCK FAST SET, and NEW GC FUJIROCK IMP available in 50 pound bulk packaging.

Characteristics and benefits:

- Bulk packaging 50 pound
- Easy to mix and flows well, yet runs less from the impression when pouring.
- Exhibits superior handling characteristics
- Increased compressive strength results in no fractures during removal of impression.
- Low setting expansion provides accurate models with smooth surfaces.
- NEW GC FUJIROCK 0.08% expansion
- NEW GC FUJIROCK FAST SET 0.08% expansion
- NEW GC FUJIROCK IMP 0.06% expansion



4+1 FREE


50 lb. Bulk Adds Great Value!

\$116⁷⁵

315020 50lb.
NEW FUJIROCK
Fast Golden Brown

Inspired by SIMPLICITY AND AFFORDABILITY

The CST® technique makes the fabrication process both fast and simple. A CST® framework can be made in approximately 30 minutes and requires no casting, scanning, or milling making it accessible to any dental professional with a light curing unit. All the usual protocols for fixed-hybrid dentures are followed which means familiar routines and processes are maintained. Not only are CST® frameworks strong and easy to fabricate, they're affordable – making a fixed denture solution more accessible to more people.



Buy a CST Intro Kit and Receive a Package of 2mm Braided Rope White or Pink!
\$79 VALUE!



CST® complete kit (CST50): \$475

Contents: 2 - 17.75" (450mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:6), 2 - 11.75" (300mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:4), 1 - 3ml syringe of FIBER FORCE® CST® Flow Pink Resin with tips, 1 - 5ml bottle of FIBER FORCE® CST® Bond Resin, 4 - FIBER FORCE® CST® Temporary Support Pillars & Instructions.

CST® refill kit (CST60): \$175

Contents: 1 - 17.75" (450mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:6), 1 - 11.75" (300mm) FIBER FORCE® CST® Pink Fiber (Hybrid 1:4), 2 - FIBER FORCE® CST® Temporary Support Pillars & Instructions.

A company on the move!



Founded in 2007 in Montreal Canada, Dental Wings Inc. is today a world leader in dental CAD/CAM technology. From its early days the company became known for its state of the art dental model scanners and comprehensive CAD design software, DWOS. Today DWOS is used by more than 7000 customers worldwide who together produces millions of digital prosthesis per year. The company has experienced sustained growth and now counts 85 employees with offices in Montreal Canada, Lyon France and Berlin Germany.

Straumann and 3M adopt DWOS Platform

Announced at the IDS 2011, the partnerships with Straumann and 3M ESPE was a watershed moment for Dental Wings. With two of the dental industry's leading companies adopting its DWOS software platform, Dental Wings reputation for meeting the highest standards of quality and precision in CAD-CAM was secured. The rigorous and demanding process of meeting the requirements of both Straumann and 3M ESPE has afforded Dental Wings a robust technological backbone that bodes well for its future.



Open, Connected and Trusted

Through collaborative development Dental Wings has established powerful partnerships with leading dental technology and material manufactures to create validated workflows using the open STL standard. Labs who own a Dental Wings system automatically have access to materials such as 3M and VITA as well as integrated workflows such as Straumann CARES, 3M Lava, Atlantis and BEGO. Importantly, the DWOS platform allows for the import of STL files from intra-oral scanning devices such as the 3M True Definition, iTero and E4D.



Integrate Digital Dentistry

With the acquisition of coDiagnostiX from Straumann in January 2013, DWOS has become the first fully open and integrated dental software suite that enables powerful cross-application integration. Now the clinician-lab team can plan prosthetic-driven guided implant surgery, design and fabricate surgical drill guides as well as the provisional restorations, all on the same open platform in one design session. Dental Wings received FDA approval in July of 2013 for its newly released coDiagnostiX 9 that enables laboratories to economically produce surgical drill guides locally or through service centers on Stratasys Objet 3D printers.

Dental Wings is a company on the move, committed to creating innovative digital dental solutions that deliver value and peace of mind to its user community.



7SERIES



3SERIES



Price Lease from **\$630⁰⁰**
per month
\$29,500

Price Lease from **\$405⁰⁰**
per month
\$19,000

Available software applications



guided surgery



implants



crown & bridge



model builder



full dentures



partial frameworks



bite splints



orthodontic archiving



SPECIAL OFFER

Receive 2 extra years of free support & upgrade after first year of warranty!

\$5000 value
7Series

\$3000 value
3Series



Unique Brush - The Truly Esthetic Brush Ceramic Processing "par excellence"

With newly developed matt black bristles and unique design of each brush - even more benefits when processing ceramic materials.

- Perfectly shaped brush tips enables accurate placement of the precious ceramic materials.
- Superior retention of moisture for prolonged modeling.
- High elasticity of the brush hair so that large quantities of ceramic material can be picked up.
- Spatula function for layering small gaps.
- Individual color of each brush renders the tool unique.
- Ergonomic shape enables fatigue-free working.

Special
\$69⁹⁵ Set (4)
Reg. \$87.89



breident

Product name	Size	VPE	REF
Unique Brush Set	1, 4, 6, 8	1 piece each	390 USET 1
Unique Brush	1	1 piece	390 U001 0
Unique Brush	4	1 piece	390 U004 0
Unique Brush	6	1 piece	390 U006 0
Unique Brush	8	1 piece	390 U008 0
Unique Brush	BigBrush	1 piece	390 U008 B

Breident Exacto-Form: Polyethylene Modelling Resin

The perfect base for precision crown and bridge construction. Exakto-Form's minimal shrinkage ensures a high accuracy of fit. The excellent flowability avoids bubbles and forms accurate margins.

Available in five colours:

breident

Isolating liquid CODE: 52000210
\$31.30 / 125 ml



- green (52000280)
- blue (52000281)
- ivory (52000282)
- yellow (52000284)
- grey (52000283)

Special
\$99⁹⁵
Reg. \$110/box

kit contains 6 tins each of component A & B

Acrylic polishing set

350 0084 0

- 1 x 150 g Abraso-Star K50 slightly abrasive
- 1 x 500 g Pumice polishing rpaste
- 1 piece Abraso-Soft Acryl
- 1 piece Abraso-Buff Acryl
- 1 piece High luster buff Acryl

Price
\$85⁶⁵

breident



The pumice polishing paste diffuses into the brushes and allows particularly long prepolishing.



The fine abrasive components of the pumice polishing paste simplify careful polishing of the acrylic denture material.

Abraso-Soft Acryl
Ø 80 mm
1 piece

Abraso-Buff Acryl
Ø 80 mm
1 piece

High luster buff Acryl
Ø 100 mm
1 piece

Pumice polishing paste
for polishing of acrylics
and metal
3 x 500 g

Titanium Finishing Set

NEW!

breident

Finish titanium rationally using cutters, polishers, brushes and pastes developed specifically for use on titanium.

Assortment:

- 1 Abraso-Soft Metal csf 2/80
- 1 Abraso-Star asg, 50 ml 1 cotton buff, for handpiece
- 1 beech wood stand 8Bo/HP
- 1 Diatit tungsten carbide cutter D001 KT 14
- 1 Diatit tungsten carbide cutter D194 KT 40
- 1 Diatit tungsten carbide cutter D194 KT 50
- 1 Diatit tungsten carbide cutter 198 KT 23
- 1 high-lustre buff, metal, 50/100
- 1 round brush zwm db 19 Ø
- 1 Titapol pre-polishing paste, 150g.
- 1 Titapol pre-polishing wheel & 1 Titapol pre-polishing cylinder.



35000890

Price
\$317⁵³

Additional VKS-OC-SG Mini Assortment Kits

Need to order our popular VKS System components in smaller quantities? These convenient assortment kits are now available in 2.2 mm and 1.7 mm!



43005350



43006720



43005520



43006780



43005360



43006630

breident

B 43005350 VKS-OC-UNI 2.2 mm **\$76.65**

10pc. Kit includes 2 each of green, yellow, red OC matrices, 2 blackout discs and 2 universal males.

B 43005520 Root Caps Kit 2.2 mm **\$129**

12pc. Kit contains 2 each of green, yellow, red OC matrices, 2 blackout discs, 2 universal males and 2 metal housings.

B 43055360 VKS-SG 2.2 mm **\$68**

8 pc. Kit includes 2 each of green, yellow, red SG matrices, 2 SG males.

B 43006720 VKS-OC-UNI 1.7 mm **\$78.15**

10pc. Kit includes 2 each of green, yellow, red, 2 blackout discs and 2 uni males.

B 43006780 VKS-OC 1.7 mm **\$126.75**

12 pc. Kit includes 2 each of green, yellow, red matrices, 2 blackout discs, 2 universal males and 2 metal housings.

B 43056630 VKS-SG 1.7 mm **\$68**

8 pc. Kit includes 2 each of green, yellow, red SG matrices, 2 SG males.

Breident Diagen Turbo Set

Diagen turbos are very aggressive sintered diamonds that work effectively on zirconium oxides such as Procera™, Emax™ and all other ceramic materials. Each stone is molded to provide optimum grinding properties and cuts effectively down to the last layer.


34000200
34000205

Diagen Turbo Grinder Set

breident

Special
\$199⁰⁰
Reg. \$225

34000200



34000205

Special
\$239⁰⁰
Reg. \$275

STOCK UP FOR THE FALL - Everyday Best Sellers

Molloplast B 170g.



#03001

Reg. \$249.00

\$224.00

SAVE
\$25

Lucitone 199 25lb.



Reg. \$918.00

\$818.00

SAVE
\$100

Edenta



5+1 Burs

C/D #2 Baseplate Wax



SAVE
\$5

Reg. \$41.50

#0252

\$36.50

Elite HD (variety)



4+1 FREE

Zhermack Alginate



10 + 2 FREE

Micro-Torch



#04020

\$27.50

Denture Baths



Imprinting Now Available!

Special
.70¢

With Imprint
\$1.45

Zetalabor



3 + 1 Available in 2.6kg
5kg & 10kg

EXA-Technique



6 + 1 FREE

Cotton Buffs 4 x 42



#1180170

\$56.20 (Pkg.12)
Purchase 1 Dozen &
Receive 2 FREE Buffs!

B12 Brushes



#B12

\$15.95 (Pkg.12)
Purchase 1 Dozen &
Receive 2 FREE Brushes!

C/D Repair Liquid



Reg. \$41.95 #00323

\$35.95

Al-Cote



Reg. \$104.97 #652502

\$93.00

DENSEP Separator



SAVE
\$10

Reg. \$53.95 #555-8

\$43.95

Denture Brushes



#00002

50 Pkg.

Reg. \$65.95

\$56.95

SAVE
\$10

AcryPlus **evo**

The Newest **eVO**lution in Teeth from Ruthinium Group the Manufacturers of the Popular Acrylux Brand.

eVOlution in Aesthetics

- Cutting the neck, dentin and multilayer stratified enamel to give a natural pearly finish.
- Texturing and characterisation using cutting edge CAD-CAM technology give a natural finish to the teeth.
- Wide neck forms enable easier aesthetic adaptation of the prosthesis to the structure beneath.
- Large forms of incisors and molars to give natural beauty to the prosthesis.

eVOlution in Technology

- Resin elasticity: teeth in PMMA resin are extremely elastic and have better resistance to the mastication force.
- High abrasion resistance: the use of high molecular weight resins guarantees long-lasting resistance for the prosthesis.
- The use of PMMA as the base guarantees a perfect bond between the tooth and resin, eliminating the problem of the single teeth detaching from the prosthesis.



\$123.⁰⁰
VALUE!

Buy 2 x (1x28)
Receive 1x28 FREE!

DESIGNED TO GET THE BEST RESULTS WITH:

- COMPLETE PROSTHESES
- OVER DENTURES
- ALL ON FOUR
- BAR PROSTHESES



UPCOMING COURSES & EVENTS



Impression Taking & Soft Tissue Model Course

Hands on course using live patient and verification jig
Sept 12, 2014 - Lethbridge AB

Implant Supported Fixed Bar Prosthesis

Hands-on Lab and Clinical based training
with Live Patient experience
Sept 13 -14, 2014 - Lethbridge AB

2014 Perfecting Your Practice

The Denturist Association of Ontario PYP Conference
Sept 18 - 20, 2014 - Mississauga ON

Spectrum Vancouver

Sept 19, 2014 - Vancouver BC

Dental Technology Academy

DTA Lecture Series
Sept 20, 2014 - Scarborough ON

Expodent

Sept 25 - 27, 2014
Montebello Quebec

Denturism Implant Study Club

DISC Lecture Series
Sept 27, 2014 - Scarborough ON

Denttechnica du Quebec

Oct 8, 2014 - Montreal Quebec

Spectrum Toronto

Oct 24, 2014 - Toronto ON

Dental Technology Academy

Evaluating Investments Opportunities for Your Lab
DTA Lecture Series
Nov 08, 2014 - Scarborough ON

Swiss Denture Concept 1 - Master Dentures

Nov 10 - 11, 2014 - Mississauga On

Swiss Denture Concept 2 - Clinical Course

Nov 14 - 16, 2014 - Concord On