

STERILIZATION STATIONS



STR-2

Regular \$6,395.00
\$5,595.00
Special



STR-3

Regular \$7,195.00
\$6,395.00
Special

- Produced of CNC Manufactured .090mm DKP STEEL.
- Electrostatic powder coating paint process.
- Sections for autoclave
- Soft closing system for drawers
- Vertical opening upper cabinet with built-in dispensers

TUTTNAUER



EZ9 Automatic Autoclave

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- Highly efficiency air pump circulates air in chamber during drying cycle
- 0.2um Hepa filter (replaceable)
- Program parameters can be customized and stored
- Optional printer documents date, time, temperature and pressure
- Exceeds regulatory requirements for infection control
- Chamber dimensions: 9" x 18" (5gal)
- Overall dimensions: 21.5" x 20" x 14.4"
- Cycle time: 27min (Cold) & 13min (Hot)

EZ9 without Printer

\$5542 (#EZ9)

EZ9 with Printer

\$6294 (#EZ9P)



SPECIAL!

Buy an EZ9 Autoclave...
Get a FREE Ultrasonic Cleaner!

RENFERT

Renfert's new Automatic Denture Cleaner!

SYMPRO

High-performance cleaning device for prosthetic dentures, orthodontic appliances and splints



For additional details go to www.renfer.com

(#6500-0000)

Special

\$895.00
Reg. \$995

+ 2L Liquid FREE!

PRODUCTS BEYOND INNOVATION

VPS impression material

elite[®] HD+
HYDROPHILIC



VPS impression material that combines the effects of Nano-Technology and hydrophilia.

- High flowability due to nano-technology
- Super ability to displace moisture resulting in excellent detail reproduction
- Excellent tear strength
- Available in FAST set and REGULAR set
- Wash, Heavy Body, monophasic & Putty materials

offer 167

BUY 4 Elite[®] HD 2x 50 ml Packs **GET 1** of the same **FREE!**
or **BUY 4** Elite[®] HD Putty's **GET 1** of the same **FREE!**

offer 168

BUY 3 Elite[®] HD Maxi 6x 380 ml cartridges **GET 1** of the same **FREE!**



ALGHAMIX II allows the mechanical mixing of *alginates, plasters and silicones* for dental use

Buy 60 (1 lb) bags of any Zhermack Alginate GET an Alghamix 2 FREE!



Modern design

alghamix II

APPLICATIONS

■ Alginate mixing



■ Silicone mixing



■ Dental stone mixing



EVA Digital X-Ray System



EVA digital x-ray system delivers big value.
Backed by the proven reliability of tens of thousands of installations ranging from single chair clinics to dental chains, EVA offers a rugged sensor design coupled with the highest standard of technology.

- Exceptional Image Quality
- Ideal For Any Practice
- Full Integration
- Best Supported Sensor
- Great Value

EVA Digital X-ray Complete System

\$7695.00

Contains: EVA Sensor, EVAsoft Imaging Software, USB Station, USB Cable, Box of Sheaths, Rinn Positioning Aids, Drivers
Warranty: 2 Year Limited.



COLTENE

ParaCore

Special
\$210.15

3 Indications – 1 Material

ParaCore is a dual-cured, glass-reinforced composite resin, radiopaque composite system available in three shades. The well-known 3-in-1 material simplifies the post & core restorative technique, because it is ideal for post cementation, core-build ups and crown & bridge cementation. The use of ParaCore as a fixing and abutment material establishes an optimum "monoblock bond interface" between dentin, post and crown to produce a restoration with an unusually long life and very high strength.



#C5885



ParaPost

ParaPost® Fiber Lux / Taper Lux

Parallel and Tapered Translucent Fiber Post Systems. An alternative to metal posts when an esthetic metal-free restoration is desired

#PF170

Special
\$226.80

OneCoat 7.0

- Single component adhesive
- Light-cured self-etching
- Designed to etch, prime and bond



#C7047

Special
\$110.00



GC America's Quarterly Specials



EQUIA™ & "RESTORATIVE RENAISSANCE" SUPER DEAL

- Buy 2 EQUIA™ Fil Refills, **Get 1 FREE!*** promo 1670
- Buy \$2,500 Retail of EQUIA™, KALORE™, G-aenial™ Flo, Universal Flo or Bond, **Get 1 Kindle Fire† FREE!*** promo 1671
Retail Value \$200 via Amazon with free shipping and no tax
- Save over 30%** on EQUIA™ Introductory Packages and Starter Packages!
Contact your GC America Representative to find out how!



RESTORATIVES

Composite & Bonding

- Buy 2 KALORE™, G-aenial™ Flo, Universal Flo or Bond Refills, **Get 1 FREE!*** promo 1676
- Buy 3 GRADIA® DIRECT or GRADIA® DIRECT X Refills, **Get 1 FREE!*** promo 1677
- Buy 3 G-BOND™ Refills, **Get 1 FREE!*** promo 1678
Same type

Glass Ionomer

- Buy 3 Glass Ionomer Restorative Products, **Get a \$75 Best Buy, Target or Amazon Gift Card FREE!*** promo 1679
Excludes GC Fuji VARNISH™, EQUIA™ Coat, GC CAVITY CONDITIONER & GC SELF CONDITIONER



CEMENTS

Temporary

- Buy 2 GC Fuji TEMP LT™ Refills, **Get 1 FREE!*** promo 1673
- Buy 3 GC TEMP ADVANTAGE® Refills, **Get 1 FREE!*** promo 1672
Same size

Resin

- Buy 3 G-CEM™ Refills, **Get 1 FREE!*** promo 1674

Glass Ionomer

- Buy 3 Glass Ionomer Cement Products, **Get a \$75 Best Buy, Target or Amazon Gift Card FREE!*** promo 1675
Excludes powder/liquid products



PREVENTIVE

MI Varnish™ will be available February 2012

- Buy 2 **NEW** MI Varnish™ Packages, **Get 1 NEW GC Tri Plaque ID Gel™ Package FREE!*** promo 1669
- Buy 2 **NEW** MI Varnish™ Packages, **Get a \$25 Best Buy, Target or Amazon Gift Card FREE!*** promo 1666
- Buy 3 MI Paste™ or MI Paste Plus™ 10-Tube Packages, **Get 1 NEW MI Varnish™ Package FREE!*** promo 1668
- Buy 3 MI Paste™ or MI Paste Plus™ 10-Tube Packages, **Get 1, PLUS 1 NEW Patient Indication Card Kit FREE!*** promo 1682
- Buy 3 **NEW** GC Tri Plaque ID Gel™ Packages, **Get 1 FREE!*** promo 1683



IMPRESSION MATERIAL

- Buy 4 EXA Products, **Get 1 FREE!*** promo 1680
Excludes 32 packages
- Buy 3 Tray Products, **Get 1 FREE!*** promo 1681



EQUIPMENT

- Buy 6 Capsule Intro Kits or Refills, **Get 1 CAPSULE MIXER CM-II™ FREE!*** promo 1684
- Buy \$1,200 Light Cured Restorative Products, **Get 1 Silverlight™ FREE!*** promo 1687



MISCELLANEOUS

- Buy 3 EPITEX® Packages, **Get 1 FREE!*** promo 1685
- Buy 3 Asepsis Products, **Get 1 FREE!*** promo 1686



CENTRALIZE

Your CAD/CAM Needs!



Implantology Is In Our Future - Is It In Yours?

THE FUTURE IS OPEN CAD/CAM, LASERDENTA IS BRINGING IT TO YOU!



Laserdenta

OPENSCAN - OPENCAD - OPENMILL

Laserdenta, Inc.
Toll Free: 866-328-5011
<http://www.laserdentaincusa.com>

graphic design by Nabuki Digital Designs - www.mynabuki.com

Laserdenta's Authorized Dealer:

C/D *Central Dental*
The laboratory and denturist supply specialist.

Toll free: 1-800-268-4442

Plastic Cups



#30-5148

5 oz Cup
(Package of 1000)

Special
\$26.95
Reg. \$31.95

Headrest Covers



Special #24-1014
\$24.95
Reg. \$28.95

Special #24-1017
\$28.95
Reg. \$32.50

#24-1014 (Headrest Cover Tissue/Poly -10"x10")
#24-1017 (Headrest Cover Tissue/Poly -10"x13")

Dental Bibs

- #23-1210 - Blue
- #23-1211 - Green
- #23-1212 - White
- #23-1213 - Lavender
- #23-1215 - Yellow
- #23-1216 - Rose
- #23-1217 - Peach



Special
\$19.95
Reg. \$23.95

Sterilization Pouches



Special #P090135
\$8.95
Reg. \$10.30

Special #P140255
\$14.95
Reg. \$16.80

#P090135 - 3.25" x 5.25" Package of 200
#P140255 - 5.5" x 10" Package of 200

NEW!

Marly Skin Guard™

When rubbed on the skin forms a two-dimensional skin barrier which dries in 3 minutes. Marly Skin Guard™ dispenses as a foam and protects the skin, for at least 4 hours, against abuse from the elements by inhibiting absorption of:

- alkaline & acidic solutions.
- organic solvents & dust.
- aggressive cleansers & disinfectants.
- chemicals & lubricants.
- bodily secretions, blood, urine & feces.
- paints, inks & dyes.



\$34.95
/25¢ per Application

Micro-Kleen3

is ready to use, broad spectrum surface & equipment disinfectant cleaner. Micro-Kleen3's fast "Accelerated Kill Time Formula™", will effectively clean and disinfect surfaces such as: Dental equipment, operating lights/tables, impression materials prosthetics, counter tops, sinks, reception areas and all "HIGH CONTACT, HIGH TOUCH" environmental healthcare surfaces.



(#MKS12-024-1) 24oz. Sprayer bottle **\$9.50**
(#MKW06-100) Case of 6-Wipes (100)/Roll **\$77.70**
(#MKP04-128-1) 3.8L bottle **\$24.95**

RONCO'S PREMIUM LINES OF EXAMINATION GLOVES

NITECH®

RONCO NITECH®

- Patented blend of premium quality synthetic materials
- Powder free



Price:
\$7.95/box
\$7.55/box (case)



NITRILE

RONCO NITRILE

- Diamond textured for superior grip
- Powder free



Price:
\$11.95/box
\$10.95/box (case)



RONCO TOUCH

- Premium quality nitrile in pink colour
- Powder free



LATEX

RONCO SILKTEX

- Superior barrier protection with excellent elasticity
- Available as powder free and lightly powdered

Price:
\$8.95/box
\$7.95/box (case)
(Powder Free)

Price:
\$7.95/box
\$6.95/box (case)
(Lightly Powdered)



www.ronco.ca

RONCO
focus on quality

BORDER LOCK®

MOULDBLE DISPOSABLE IMPRESSION TRAYS

Utilizing the design concepts of the original Border-Lock® Impression Tray System, new Mouldable Border-Lock Disposable Impression Trays take impression accuracy one step further. Offering the ability to modify the trays chairside by hand, to perfectly adapt to your patient's arch shape and size, Mouldable Border-Lock trays provide a true, anatomically correct fit creating optimal dynamic pressure, which greatly reduces the risk of impression distortion. To use, simply select the appropriate tray size using the calipers and sizing chart provided. Immerse the selected tray in a water bath of 70°C for at least one minute; once softened, the tray can be adapted by hand. After adaption, the tray will assume its normal rigidity. Mouldable Border-Lock trays outperform stock metal trays and allow you to achieve results that rival those of custom trays, yet more conveniently and more efficiently. For precision impressions, this is the system to use!

Refills 5/pkg \$21.95

Edentulous Sizes: 21-25 & 31-35
Dentulous Sizes: 42-49 & 51-58

Special

Buy 5, Get 1 FREE



RUTHINIUM

Dust-Free Alginate

- Fast mixing (30 sec.)
- Setting time: 2' (with water at 23°C/73°F)
- Flavour: strawberry
- Colour: pink

Special Reg. \$9.95
\$8.95
/bag

when buying a case (12+ bags)!



#2060001



Buy 2 Bags of Alginate,
Get a Mixing Bowl
& Spatula FREE!