



CENTRAL DENTAL - WORK BENCH CATALOGUE



Tel: 1(800)268-4442 Fax: (416)694-1071 www.CentralDentalLtd.com

3420 Pharmacy Ave. Unit 3, Scarborough, ON, M1W 2P7

Technical Properties of Laboratory Workbenches

- Made of 1.2mm DKP Steel Sheet
- Metal surfaces are plated with Ferrum-phosphate (Corrosion Resistant)
- Painted by electrostatic powder paint method
- A comprehensive range of RAL colors is offered
- Work-top is made of 20mm solid compact, a synthetic material which is heat and scratch resistant
- Ergonomic lights
- Air nozzle
- Counter balance air nozzle return system
- Vacuum system built-in
- Variable speed control for vacuum system
- Wooden arm rest and protective stainless steel work sheet
- 2 slim drawers underneath the work area
- Minimum 4 drawers, depending on bench model
- Separate light switch with neon indicator
- Each work position fitted with top power receptacles
- Electric power receptacle underneath each work position
- Central Dental work benches are certified by "Electrical Safety Authority"

Optional Extras

- 25mm parquet designed hard wood worktop (beech) satin polished
- Adjustable lighting system
- Chrome legs
- Soft touch drawer system
- Wood suction attachment
- Chrome Handles
- Rounded stainless steel work plate

Terms and Conditions

- 40% Deposit with order
- 40% Upon Delivery
- 20% 30 Days from receipt



M-82 Single Laboratory Workbench

Technical properties

- Dimensions: 111.5cm x 60cm x 86cm(h)
- Weight: 85 kg
- 1250 Watt power consumption

Items Included

- 2 slim drawers, 3 medium drawers, 1 large drawer
- Full suction system with the bench attachment
- 1 air nozzle
- 2 arm rests
- Full single span fluorescent light with built in work shelf



M-82 - A Single Workbench

Technical properties

- Dimensions: 150cm x 60cm x 86cm(h)
- Weight: 85 kg
- 1250 Watt power consumption

Items Included

- 2 slim drawers, 3 medium drawers, 1 large drawer
- Full suction system with the bench attachment
- 1 air nozzle
- 2 arm rests
- Full single span fluorescent light with built in work shelf



M-82 - C Single Workbench

Technical properties

- Dimensions: 178cm x 178cm x 60cm x 86cm(h)
- Weight: 195 kg
- 1300 Watt power consumption

Items Included

- 4 slim drawers, 6 medium drawers, 2 large drawers
- 2 air nozzles
- Full suction system with the bench attachment
- 2 full span fluorescent light with built in work shelf
- 2 arm rests



M-93 Double Laboratory Workbench

Technical properties

- Dimensions: 175cm x 60cm x 86cm(h)
- Weight: 130 kg
- 2500 Watt power consumption

Items Included

- 4 slim drawers, 3 medium drawers, 1 large drawer
- 2 Full suction system with the bench attachment
- 2 air nozzle
- 4 arm rests
- 1 double span fluorescent light with built in shelf



M-96 Double Laboratory Workbench

Technical properties

- Dimensions: 215cm x 60cm x 86cm(h)
- Weight: 155 kg
- 2500 Watt power consumption

Items Included

- 4 slim drawers, 6 medium drawers, 2 large drawer
- 2 full suction systems with bench attachment
- 2 air nozzles
- 4 arm rests
- One triple span fluorescent light with built in shelf



M-100 Three Grouped Asterisk Table

Technical properties

- Dimensions: 240cm x 280cm x 86cm(h)
- Weight: 430 kg
- 1300 Watt power consumption

Items Included

- 12 slim drawers, 15 medium drawers, 6 large drawers
- Includes 46cm x 49cm cupboard 10cm below the working area
- 6 air nozzles
- 9 Electric plug sockets on working area
- 3 suction systems & 6 arm rests
- Specially designed section for a porcelain furnace & vacuum pump
- 3 way lighting system



M-SAM

Technical properties

- Dimensions: 350cm x 350cm x 86cm(h)
- Weight: 600 kg
- 8000 Watt power consumption

Items Included

- 16 slim drawers, 26 medium drawers, 4 large drawers
- Includes 46cm x 49cm cupboard 10cm below the working area
- 8 air nozzles
- 4 suction systems
- 12 Electric plug sockets on working area
- Specially designed section for a porcelain furnace & vacuum pump
- 4 arm rests
- 4 way lighting system



M-200

Technical properties

- Dimensions: 260cm x 310cm x 86cm(h)
- Weight: 550 kg
- 8000 Watt power consumption

Items Included

- 8 slim drawers, 16 medium drawers
- 4 suction system with bench attachment
- Includes 46cm x 49cm cupboard 10cm below the working area



E Series - Cabinets

Technical properties

- Can be used in conjunction with any workbench (same height)
- Made of 1.20mm DKF Steel Sheet
- Metal surfaces are plated with Ferrum-Phosphate (corrosion resistant)
- Painted by Electrostatic powder paint method & kiln dries @ 200°C
- A comprehensive range of RAL colours is offered



E-35

Technical properties

- Dimensions: 60cm x 60cm x 86cm(h)
- 2 doors and 2 shelves
- Various top choice



E-36

Technical properties

- Dimensions: 50cm x 60cm x 86cm(h)
- 1 door and 2 shelves
- Various top choice



E-39

Technical properties

- Dimensions: 50cm x 60cm x 86cm(h)
- 4 drawers
- Various top choice



E-50

Technical properties

- Dimensions: 50cm x 60cm x 86cm(h)
- 1 door and 1 shelf
- Various top choice

LAB 21 Series Plaster Solutions

Technical properties

- Made of 1.20mm DKF Steel Sheet
- Metal surfaces are plated with Ferrum-Phosphate (corrosion resistant)
- Painted by Electrostatic powder paint method & kiln dries @ 200°C
- A comprehensive range of colours is offered
- Worktop is made of 2mm stainless steel
- Units can be fitted with plaster waste chute and bin
- All units are modular and can be combined to create the optimal working area



LAB-21 A

Technical properties

- Dimensions: 50cm x 55cm x93cm(h)
- 1 door and 2 shelves
- Stainless steel worktop



LAB-21 AA

Technical properties

- Dimensions: 100cm x 55cm x93cm(h)
- 2 doors and 2 shelves
- Stainless steel worktop



LAB-21 AD

Technical properties

- Dimensions: 100cm x 55cm x93cm(h)
- 1 door, 2 shelves and 5 drawers
- Stainless steel worktop



LAB-21 AK

Technical properties

- Dimensions: 100cm x 55cm x93cm(h)
- 1 door and 2 shelves
- Corner element
- Stainless steel worktop



LAB-21 B

- Dimensions: 50cm x 55cm x93cm(h)
- 1 door, 1 shelf and 1 drawer
- Stainless steel worktop



LAB-21 BB

- Dimensions: 100cm x 55cm x93cm(h)
- 2 doors, 1 shelf and 2 drawers
- Stainless steel worktop



LAB-21 C

- Dimensions: 50cm x 55cm x93cm(h)
- 3 plaster drawers and plaster basins
- Stainless steel worktop



LAB-21 D

- Dimensions: 50cm x 55cm x93cm(h)
- 5 drawers
- Stainless steel worktop



LAB-21 DA

- Dimensions: 100cm x 55cm x93cm(h)
- 1 door, 2 shelves and 5 drawers
- Stainless steel worktop



LAB-21 DD

- Dimensions: 100cm x 55cm x93cm(h)
- 10 drawers
- Stainless steel worktop



LAB-21 E1

- Dimensions: 100cm x 55cm x 93cm(h)
- 2 doors and plastic trap
- 35x40 x18cm(h) small size sink
- Mechanic air system and air syringe
- Stainless steel worktop



LAB-21 E2

- Dimensions: 100cm x 55cm x 93cm(h)
- 2 doors and plastic trap
- 35x70 x18cm(h) small size sink
- Mechanic air system and air syringe
- Stainless steel grate for model trimmer
- Stainless steel worktop



LAB-21 G1

- Dimensions: 50cm x 55cm x 93cm(h)
- 1 door and 1 shelf
- Waste hole and waste bin
- Stainless steel worktop



LAB-21 GA

- Dimensions: 100cm x 55cm x 93cm(h)
- 1 door and 1 drawer
- Plaster storage cabinet
- Waste hole and waste bin
- Stainless steel worktop



LAB-21 GB

- Dimensions: 100cm x 55cm x 93cm(h)
- 1 door and 1 shelf
- 2 plaster drawers
- Waste hole and waste bin
- Stainless steel worktop



LAB-21 GC

- Dimensions: 100cm x 55cm x 93cm(h)
- 1 door and 1 shelf
- 3 plaster drawers and plaster basin
- Waste hole and waste bin
- Stainless steel worktop



LAB-21 AS1

- Dimensions: 50cm x 32cm x 60cm(h)
- 1 door and 1 drawer
- Wall cabinet



LAB-21 AS2

- Dimensions: 100cm x 32cm x 60cm(h)
- 2 doors and 1 shelf
- Wall cabinet



Accessories




ER-L 02 Ceramist Light

- 2x18 & 2x36 Wall
- Anti-Dazzle aluminum grid



ER-L 03

- Adjustable day light
- 40 Watt & 3000 lm



**ER-0013.W
Wood Attachment**



**ER-0014
Air Syringe**



**ER-0015
Retractable
Mechanism**



**ER-0016
Arm Support**



**ER-0013
Plastic Attachment**



Orthopedic Chair

- Adjustable at 3 positions



RAL Colour Options:

Colour Options:



These should be treated as an approximate guide as all printers and screens will replicate these colours differently

RAL1000	RAL1001	RAL1002	RAL1003	RAL1004	RAL1005	RAL1006	RAL1007	RAL1011	RAL1012
RAL1013	RAL1014	RAL1015	RAL1016	RAL1017	RAL1018	RAL1019	RAL1020	RAL1021	RAL1023
RAL1024	RAL1027	RAL1028	RAL1032	RAL1033	RAL1034	RAL2000	RAL2001	RAL2002	RAL2003
RAL2004	RAL2008	RAL2009	RAL2010	RAL2011	RAL2012	RAL3000	RAL3001	RAL3002	RAL3003
RAL3004	RAL3005	RAL3007	RAL3009	RAL3011	RAL3012	RAL3013	RAL3014	RAL3015	RAL3016
RAL3017	RAL3018	RAL3020	RAL3022	RAL3027	RAL3031	RAL4001	RAL4002	RAL4003	RAL4004
RAL4005	RAL4006	RAL4007	RAL4008	RAL4009	RAL4010	RAL5000	RAL5001	RAL5002	RAL5003
RAL5004	RAL5005	RAL5007	RAL5008	RAL5009	RAL5010	RAL5011	RAL5012	RAL5013	RAL5014
RAL5015	RAL5017	RAL5018	RAL5019	RAL5020	RAL5021	RAL5022	RAL5023	RAL5024	RAL6000
RAL6001	RAL6002	RAL6003	RAL6004	RAL6005	RAL6006	RAL6007	RAL6008	RAL6009	RAL6010
RAL6011	RAL6012	RAL6013	RAL6014	RAL6015	RAL6016	RAL6017	RAL6018	RAL6019	RAL6020
RAL6021	RAL6022	RAL6024	RAL6025	RAL6026	RAL6027	RAL6028	RAL6029	RAL6032	RAL6033
RAL6034	RAL7000	RAL7001	RAL7002	RAL7003	RAL7004	RAL7005	RAL7006	RAL7008	RAL7009
RAL7010	RAL7011	RAL7012	RAL7013	RAL7015	RAL7016	RAL7021	RAL7022	RAL7023	RAL7024
RAL7026	RAL7030	RAL7031	RAL7032	RAL7033	RAL7034	RAL7035	RAL7036	RAL7037	RAL7038
RAL7039	RAL7040	RAL7042	RAL7043	RAL7044	RAL7045	RAL7046	RAL7047	RAL8000	RAL8001
RAL8002	RAL8003	RAL8004	RAL8007	RAL8008	RAL8011	RAL8012	RAL8014	RAL8015	RAL8016
RAL8017	RAL8019	RAL8022	RAL8023	RAL8024	RAL8025	RAL8028	RAL9001	RAL9002	RAL9003
RAL9004	RAL9005	RAL9006	RAL9007	RAL9010	RAL9011	RAL9016	RAL9017	RAL9018	



 **Central Dental**
The laboratory and denturist supply specialist.

Tel: 1(800)268-4442 Fax: (416)694-1071 www.CentralDentalLtd.com

3420 Pharmacy Ave. Unit 3, Scarborough, ON, M1W 2P7