

Summer 2026

Flyer

Promotions in effect July 1st - September 30th 2026



Scan to place your order or speak to a representative or call toll-free:

1-800-268-4442

C/D *Central Dental*
The laboratory and denturist supply specialist.

Endless possibilities

Blue
Line®



SR Vivodent® DCL

- 42 individual anterior tooth moulds
- Made from durable, highly cross-linked DCL material
- Contemporary shade range



LEARN
MORE

Ivoclar Academy

Removable Prosthetics

Markus Fischer, RD, RDT, PID

BPS® Certification Course

Distinguish yourself in the field of Removable Prosthetics by becoming BPS Certified!

Course Description

The Bio-functional Prosthetic System (BPS) is where aesthetic brilliance meets function. BPS precision dentures are the result of a coordinated system of products, instruments, set-up techniques and training. The goal is to create the highest quality denture available. The BPS brand is available exclusively to certified BPS clinics. This two-day course consists of clinical and laboratory demonstrations and hands on processes. The course curriculum concentrates on the design and fabrication of the gnathometer M impression trays and the bio-functional set-up of an upper and lower full denture. The use of BPS clinical and technical approaches will enhance the final outcome. Participants will learn comprehensive clinical and technical training on the techniques required to create precision dentures.

Course presented in partnership with Central Dental.

When: September 24–25, 2026
8:30 am–5 pm EST

Where: Ivoclar Academy Canada
1-6600 Dixie Road
Mississauga, ON L5T 2Y2

CE: 14 CE Credits

There are two fee options for this course:

Course only: \$1399

– OR –

The Complete Bundle: \$3200

⊕ Includes a Stratos® 300
Articulator (calibrated) *Reg.* \$2,495



REGISTER
HERE



Instructor

Markus Fischer, RD, RDT, PID



Markus was born and raised in Germany and immigrated to Canada in 1986. He began his career as a dental technician in 1993 and became licensed in 1996. A year later, he started his own orthodontic dental lab and became a dental technician instructor at Vancouver Community College. Eventually he saw the light and realized he wanted to become a denturist and received his license in 2006. For the past five years his quest to make better dentures has taken him to Tokyo Japan, Leicester England, Germany, Toronto and Chicago where he studied and applied Dr. Jiro Abe's lower suction denture technique (SEMCD). He was the first Canadian denturist to be trained by Dr. Abe, is one of four Canadian SEMCD/Lower Suction Denture Instructors and provides training at his denture clinic in Mission, BC.



Ivoclar Vivadent, Inc.
Nationally Approved PACE Program Provider for RDT/RD/CDentist.
Approved does not imply endorsement.
© 2024 by I.V.I. 2024/09/15
Provider #18-1998


BEST
VALUE
SAVE
\$694

bredent

Abraso-Sil Acryl

This buff consists of a nonwoven fibre fabric in the centre between two layers of silicone-coated cotton fabric.

On the outside there are two rows of bleached Chungkink bristles. This brush absorbs polishing paste or pumice and releases it slowly for efficient prepolishing.



#3500093

Ø 80 mm
1 piece

8+1

ROPAK/Kompaktopaker UV



Ready to use Light Cure Opaker!

The ready-to-use alternative for convenient coating of CoCr objects.



Apply Ropak Kompaktopaker with the integrated brush directly onto the clean metal surface.

Reg. \$150.45
\$120.00

Ropak Kompaktopaker (Pink) 20ML 54000133

Reg. \$82.70
\$72.00

Kompaktopaker Zahnfarben (Tooth Colour) 10ML 54000105

Bredent Diagen Turbo

Diagen turbos are very aggressive sintered diamonds that work effectively on zirconium oxides such as Procera™, Emax™ and all other ceramic materials. Each stone is molded to provide optimum grinding properties and cuts effectively down to the last layer.



\$99.35

Cylinder
Ø 6.5 x 13 mm
2 pieces
REF 340 0017 0



\$70.20

Cone
Ø 3.5 x 11 mm
2 pieces
REF 340 0015 0



\$77.75

Inverted cone with recess - Ø 12 x 6 mm
1 piece
REF 340 0024 0



\$74.50

Lens
Ø 22 x 2 mm
1 piece
REF 340 0021 0




\$74.50

Disc
Ø 22 x 2 mm
1 piece
REF 340 0022 0



\$99.35

Disc
Ø 22 x 4,5 mm,
1 piece
REF 340 0019 0



6+1

of equal or lesser value

Universal Screwdriver Set

Screwdriver set for 98% of all screws available on the market. To be inserted into the torque ratchet, adjustable from 10 to 40 Ncm. This way screws can be turned in correctly and safely.

Universal screwdriver set with instruments
REF 310 0001 2

\$2195.00

Reg. \$2211.95



Torque Ratchet

Torque adjustable from 10 to 40 Ncm.

\$549.00



Torque ratchet
REF 330 0115 5



RODIN[®]
3D Resin

A Pac-Dent Brand

RODIN TRUMODEL RESIN

Precision. Stability. Confidence.
Every Print.

Rodin TruModel Resins are engineered for accuracy and reliability in crown, bridge, and diagnostic applications. Developed with FDA-listed, VOC-free chemistry, TruModel delivers consistent performance you can trust.

KEY FEATURES

Regulatory-Grade Assurance

FDA-listed formulation with VOC-free chemistry ensures verified safety, quality, and compliance for professional use.

High Dimensional Accuracy


Maintains true-to-design fits with (+/-) 50 µm on selected printers.

Matte, Scanner-Friendly Surface

Enhances margin visibility and scanning precision—no spray required.

Long Term Stability

Resists warping, chalking, and yellowing for dependable accuracy over time.



3+1
on 4 litre size

PRODUCT	ITEM #	SIZE	PRICE
Rodin TruModel Gray	#24066	1 kilo	\$89.95
Rodin TruModel Gray	#24075	4 kilo	\$320.00
Rodin TruModel Prime (Tan)	#24067	1 kilo	\$89.95
Rodin TruModel Prime (Tan)	#24076	4 kilo	\$320.00
Rodin TruModel Sand	#24065	1 kilo	\$89.95
Rodin TruModel Sand	#24074	4 kilo	\$320.00

Aidite Biomic™ Esthetic Solutions

Two separate solutions for esthetic finishing and zirconia bonding

BIOMIC™ STAIN/GLAZE



Biomic Stain/Glaze

Specially developed for CAD/CAM all-ceramic restorations. Allows dental technicians to use a fast, simple method to finish their work, achieving restorations that replicate natural teeth.

KEY FEATURES

- **CAD/CAM All-Ceramic Restorations**
Specially developed for CAD/CAM all-ceramic restorations.
- **Fast, Simple Method**
Designed to help dental technicians finish their work efficiently.
- **Replicates Natural Teeth**
Helps achieve restorations that replicate natural dentition.
- **2D Paste Colors**
Designed to match the aesthetics of natural teeth.
- **3D Paste Effects**
3D paste effects are added; combining 2D and 3D paste can maximize esthetic results.

IDEAL FOR

Natural-looking CAD/CAM all-ceramic restorations.

LiSi CONNECT™



LiSi Connect

Sprayed on the zirconia bonding surface. After sintering, LiSi Connect crystallizes into a lithium disilicate coating to complete surface modification of the zirconia.

KEY FEATURES

- **Lithium Disilicate Coating**
Crystallizes on zirconia after sintering to complete surface modification.
- **Glass Ceramic-Like Bonding**
Treated zirconia has the same clinical bonding effect as glass ceramics.
- **Tightly Combined Coating**
LiSi Connect coating is tightly combined with zirconia.
- **Ultra-Thin Micron-Level Layer**
Does not affect clinical placement.
- **Durable Adhesive Bond**
Allows long-time stable and durable adhesive bonding.
- **Compatible With Zirconia**
Can support zirconia veneer bonding and is compatible with all zirconia materials.

IDEAL FOR

Zirconia surface modification and veneer bonding.

WHY CHOOSE AIDITE BIOMIC SOLUTIONS?

NATURAL ESTHETICS

Replicates natural teeth with 2D paste colors and 3D effects.

ZIRCONIA BONDING

Lithium disilicate coating helps zirconia bonding performance.

EFFICIENT WORKFLOW

Fast, simple methods support efficient lab finishing.

Aidite Aizir™ Zirconia

Premium Multi-Layer Zirconia for Exceptional Strength & Natural Aesthetics



ONE DISC. MULTIPLE INDICATIONS. EXCEPTIONAL RESULTS.

RECOMMENDED CLINICAL APPLICATIONS



Anterior Crowns



Posterior Crowns



Anterior Bridges



Posterior Bridges




Veneers



Implant Restorations

TRANSLUCENCY & FLEXURAL STRENGTH



Natural Multi-Layer Gradient

Data based on internal testing

KEY BENEFITS

53% Incisal Translucency
Natural light transmission for lifelike esthetics

43% Cervical Translucency
Optimal opacity for shade stability and masking

800–1250 MPa Strength
High flexural strength for durable restorations

Multi-Layer Gradient
Smooth transition from cervical to incisal

Versatile Applications
Ideal for crowns, bridges, veneers and implants

Excellent Milling & Sintering
Consistent results with reliable performance



A7 Fix

P R A C T I C A L

The articulator A7 Fix brings its fixed measurements to medium standards. This allows a greater practicality and speed in the execution of the works. Unique in the world with magnetic condylar guidance (patented inlaid magnetic disc), it provides stability in joint movements, eliminating springs and elastics.

- Semi-adjustable / Arcon
- Fixed Intercondylar distance (110mm)
- Condylar guide fixed on the average: 30°
- Bennett angle fixed on the average: 15°
- Condylar Curve Guide
- Motion stabilizer system with magnet (magnetic)
- Upper Branch support pin in the open position
- Central locking

INTRODUCTORY
PRICE

\$139.95

MAGNETIC MOTION STABILIZER

Patented magnetic system for smooth, stable movements.

CENTRAL LOCKING

Quick and secure central locking system.

UPPER BRANCH SUPPORT PIN

Provides safe support in the open position.

STABLE AND RELIABLE

Engineered for durability and everyday use.

#AAR60257

A7 Plus

Revolutionary model, the A7 Plus is the semi-adjustable articulator with highest approval by the most renowned universities in the country. It is ideal for academic and professional work.

Sturdy structure, design that ensures great stability and visibility.

- Semi-adjustable / Arcon
- Fixed Intercondylar distance (110mm)
- Curved condylar guide
- Angle adjustment of the Condylar and Bennett guide
- Central Locking
- Motion stabilizer system with silicone joint
- Upper Branch Support Pin in the open position

INTRODUCTORY
PRICE

\$295.00



MAGNETIC MOTION STABILIZER

Patented magnetic system for smooth, stable movements.

CENTRAL LOCKING

Quick and secure central locking system.

UPPER BRANCH SUPPORT PIN

Provides safe support in the open position.

STABLE AND RELIABLE

Engineered for durability and everyday use.

#AAR60254

FEATURE	A7 FIX	A7 PLUS	BENEFIT
Intercondylar Distance	Fixed (110mm)	Fixed (110mm)	Standardized mounting for consistency
Condylar Guide	Fixed (30°)	Adjustable	Adjustability for greater accuracy
Bennett Angle	Fixed (15°)	Adjustable	Adjustability for greater accuracy
Stabilizer System	Magnetic	Silicone Joint	Smooth, stable condylar movements
Upper Branch Support Pin	Yes	Yes	Safe support in open position
Central Locking	Yes	Yes	Quick and secure locking
PRICE	\$139.95	\$295.00	<i>Introductory Pricing</i>

enbio.

ENBIO AUTOCLAVES

Advanced sterilization. Swiss design. Built for your practice.

ENBIO S

Class B Autoclave

Item #: 181133460US

Price: **\$5,195.00** – Special price



ENBIO PRO

NEW

Class B Autoclave

Item #: 1220000

Price: **\$7,295.00** – Special price



Unwrapped
FAST / 273°F

7 min

Wrapped
Vacuum / 273°F

15 min

Wrapped
Vacuum / 250°F

45 min

Unwrapped FAST
273°F / 134°C

9 min


Wrapped Class B
273°F / 134°C

15 min

Wrapped Class B
250°F / 121°C

45 min

CAPACITY & LOAD SPACE



Load space
11 x 7 x 1.3 in

Max. packages amount (circa) **8**
Max. recommended batch weight **1.1 lb**
Process chamber volume **2.7 l**

CAPACITY & LOAD SPACE



Load space
11.2 x 7.5 x 2.8 in

Max. packages amount (circa) **16**
Max. recommended batch weight **1.8 lb (0.8 kg)**
Process chamber volume **5.3 l**

DIMENSIONS



22 in (56 cm)



9.8 in (25 cm)

6.3 in
(16 cm)



22.6 in (57.5 cm)



10.6 in (27 cm)

7.9 in
(20 cm)

PRESSURE & TEMPERATURE

Working pressure 2.1 Bar / 30.5 psi
Max. pressure 2.45 Bar / 33.4 psi
Max. process temperature 137°C / 278.6°F

ELECTRICAL

Max. power supply 110-120V / 60Hz
Nominal power 3.25 kW
Max. power consumption 15 A

MISC.

Ingress protection IP20
Data archiving USB
Water conductivity* < 15 µS/cm

* - demineralized water only

PRESSURE & TEMPERATURE

Working pressure 2.1 Bar / 30.5 psi
Max. pressure 2.45 Bar / 33.4 psi
Max. process temperature 137°C / 278.6°F

ELECTRICAL

Max. power supply 110-120V / 60Hz
Nominal power 3.25 kW
Max. power consumption 15 A

MISC.

Ingress protection IP20
Data archiving USB
Water conductivity* < 15 µS/cm

* - demineralized water only

IMPRESSION MATERIAL PROMOTION

HUGE Perfit PVS

**3+1
FREE**



- PER003** Perfit Light Body Normal Set **PER006** Perfit Heavy Body Normal Set
PER004 Perfit Light Body Fast Set **PER008** Perfit Super Light Body Regular Set
PER005 Perfit Heavy Body Fast Set **PER009** Perfit Monophase Normal Set

Ivoclar Virtual PVS

**3+1
FREE**



HUGE Perfit Bite Registration

PER007

**4+1
FREE**



Zhermack Occlufast Rock Bite Registration

C200726



Central Dental Mixing Tips

**5+1
FREE**



PART #	DESCRIPTION	PRICE
YM110PB	Pink	\$33.28
YM120GB	Teal	\$33.25
YM100YB	Yellow	\$33.25

Buy any 5, Get 1 Free! **(5+1 Free)**

HUGE Perfit Tray Adhesive

PER010



Ensures strong bonding between VPS impression materials or plastic trays. Easy to use- just wait 3 minutes after applying before loading the impression material.

IMPRESSION TRAYS

Precision Fit. Patient Comfort. Reliable Results.

High quality trays designed for accurate impressions, easy handling, and exceptional patient comfort. Available in a variety of styles to meet your clinical needs.



KEYSTONE EDENTULOUS DISPOSABLE TRAYS

Single-use, perforated edentulous impression trays that can be easily drilled for implant impressions. These unique and economical trays allow for maximum material retention and accuracy without the use of an adhesive. High impact styrene provides excellent impressions with outstanding detail.



Perforated design for superior material retention



Easily drilled for implant impressions



High impact styrene for strength and accuracy

DESCRIPTION	PRICE
Package (12) Assorted Edentulous Perforated Trays	\$48.95



BUY 5 PACKAGES, RECEIVE 1 FREE!



GC DISPOSABLE IMPRESSION TRAYS

Spacer Impression Trays are disposable trays with a unique zigzag spacer bar, designed into the occlusal surface, which allow you to seat loaded impression trays to the correct depth, time after time. Available in either perforated or solid styles to accommodate the majority of preferred techniques.



Unique zigzag spacer bar for consistent depth



Perforated or solid styles



Disposable convenience saves time and effort

DESCRIPTION	PRICE
Package (12)	\$24.53
Package (72)	\$141.53



BUY 5 PACKAGES, RECEIVE 1 FREE!



ASA STAINLESS STEEL IMPRESSION TRAYS

High quality stainless steel trays for accurate impressions and long-lasting performance.



Only available in kit of (8)



Small, Medium and Large sizes



Upper and Lower Arch



Anatomical and ergonomic design



Autoclavable stainless steel construction

DESCRIPTION	PRICE
#S2804-8 Set of (8) Perforated Dentulous Trays	\$133.65 SAVE \$30
#S2806-8 Set of (8) Perforated Edentulous Trays	\$194.80 SAVE \$50
#S2807-8 Set of (8) Solid Dentulous Trays	\$133.65 SAVE \$30

Scheftner Premium Partial Denture Alloys

MODELSTAR S & STARBOND CO



Manufactured in Germany

Scheftner’s premium cobalt-chromium alloys provide exceptional strength, corrosion resistance, castability and finishing characteristics for removable partial dentures, clasp-retained prostheses and combined restorations.

MODELSTAR S

- Universal partial denture alloy
- Excellent for 3D printed pattern casting
- Easy finishing & polishing
- Laser weldable
- Nickel, beryllium, cadmium & lead free

STARBOND CO

- Bondable CoCr alloy
- Ideal for clasp-retained prostheses
- Combined techniques & ceramic restorations
- Spring-hard properties
- Laser weldable

Technical Comparison

Specification	ModelStar S	Starbond Co
Composition	Co 62.7% / Cr 29% / Mo 6%	Co 62% / Cr 30% / Mo 5.5%
UTS	770–780 MPa	760 MPa
Yield Strength	510–600 MPa	500 MPa
Elastic Modulus	200–230 GPa	250 GPa
Hardness	355 HV10	380 HV10
Density	8.4 g/cm ³	8.1 g/cm ³
Laser Weldable	Yes	Yes
ISO Classification	Type 5	Type 5



Flexing to meet every shade needs



PMMA Discs (Shaded and Multi-layer)

2+1

BUY 2 PMMA Discs,
GET 1 FREE.

Promo code: **CPMMAFY26**

Dentsply Sirona offer practices and laboratories a continuously growing variety of materials for economical and precise fabrication of clinically flawless and esthetically high-quality prosthetics.



Cercon® Discs

3+1

BUY 2 Cercon Discs,
GET 1 FREE.

Promo code: **DSC26-Q3-1**

20+5

BUY 20 Cercon Discs,
GET 5 FREE.

Promo code: **CCER20FY26**



Lucitone Digital Fit™ Discs

3+1

BUY 3 Lucitone Digital Fit Discs,
GET 1 FREE.

Promo code: **CLFD31FY26**

20+5

BUY 20 Lucitone Digital Fit Discs,
GET 5 FREE.

Promo code: **CLFD20FY26**

ASTRON CLEARsplint®

THERMO-FORM VACUUM FORMING MATERIAL

THE ORIGINAL

CLEARsplint®

Available in CAD/CAM
Disc or Powder & Liquid



SPECIAL!

CLEARsplint Flexible Super Bulk Kit

\$593.00

PCSD981611 – Disc 98 x 16 mm

\$99.00

PCSD9811 – Disc 98 x 20 mm

\$99.00

- ✓ Exceptional Accuracy
- ✓ Dimensional Stability
- ✓ Patient Adjustable Comfort

KEYSTONE INDUSTRIES

4+1 FREE!



BUY 4 MOUTHGUARDS – GET 1 FREE! Pkg (12)

DETAX

HIGH PERFORMANCE DENTAL MATERIALS

FREE BUR!



#03001 MOLLOPLAST B 170 g

BUY 1 MOLLOPLAST B

RECEIVE FREE TRIMMER BUR

#7275.060

DEDECO LATHE BURS



9071

9072

9073

9074

9075

9081

9082

9083

9084

9085



FREE HOLDER!

FREE FOREDOM BUR HOLDER

(ITEM# 10080)

BUY ANY 4 DEDECO LATHE BURS

AND RECEIVE A FREE

FOREDOM BUR HOLDER (ITEM# 10080)!

VITA

VITA LAB DEALER PROMOTIONS

July - September 2026

FURNACES



Buy an VITA V-60® i-Line Plus and Vacuum Pump, Get a VITA AKZENT® Plus Mini Kit FREE!
Promo Code: QFV12633

SHADE



Buy an VITA Easyshade® V or LITE, Get a VITA® Toothguide 3D-MASTER FREE!
Promo Code: QES2631

CHARACTERIZATION



Buy 3 VITA AKZENT® LC Glaze, Get 1 FREE!
Promo Code: QAKL2631




Buy 5 VITA AKZENT® LC Chroma Stains, Get 1 FREE!
Promo Code: QAKL2631



Buy an VITA AKZENT® LC Gingiva Kit, Get A Bottle of AKZENT® LC Glaze FREE!
Promo Code: QAKL2632

TEETH



Switch to **VITA Premium Teeth with a minimum \$3,500 Stock Purchase**, get a Corresponding VITA Living Mould Guide free!
Promo Code: QTHE2632

Switch to VITA MFT Teeth with a minimum \$1,500 Stock Purchase, Get a VITA Living Mould Guide and VITA TOOTHGUIDE 3-D MASTER® free!
Promo Code: QTM2632



© 2026 VITA North America. Valid on purchases made July 1 - Sept 30, 2026 only. Please allow 2 weeks for processing. Offers valid in Canada only. Offers may not be combined with any other promotion or discounts. To redeem, email invoice to promo@vitanorthamerica.com no later than Oct 15, 2026, along with requested free goods. VITA North America reserves the right to discontinue or change the offers at any time without notice. No cash value.

800-263-4778

www.vitanorthamerica.com

SHINING 3D
DENTAL



Aoralscan Elf

Smarter Occlusion

Light • Small • Powerful



Call for
a Demo
800.268.4442
\$12,600