

Fall 2025

Flyer

Promotions in effect Oct 1st - December 31st, 2025



Scan to place your order or speak to
a representative or call toll-free:

1-800-268-4442

C/D *Central Dental*
The laboratory and denturist supply specialist.

Endless possibilities

Blue
Line®



SR Vivodent® DCL

- 42 individual anterior tooth moulds
- Made from durable, highly cross-linked DCL material
- Contemporary shade range



REQUEST
FREE SAMPLE

ivoclar.com/sr-vivodent-dcl

ivoclar.com
Making People Smile

ivoclar

20%
OFF

Precision at the push of a button

IvoBase® System

Denture base material with fully
automated injection process

- High precision due to shrinkage compensation
- Fully automated process
- Predosed units for an accurate mixing ratio



LEARN
MORE



Scan to Schedule
a Virtual Demo!



Empower Your Workflow with Rodin and SOL Plus



High Production Volume Meets Precision

- Prints up to 28 full arches vertically or 8 full arches flat
- 45 µm XY pixel size creates smooth and fine prints reliably



Enhanced Print Performance

- Built-in heating rod for optimal print results
- Heats from 63°F to 83°F in 5 minutes



Validated, Versatile Solution for Dentistry

- Validated with industry-leading 3D printing resins for FDA, Canada Health and CE approved workflows



Concierge Service Included

- Remote onboarding, training and support throughout 2-year warranty



- SOL Plus with Concierge Service
- CLEANI with Dual Tanks
- CURIE Plus with N2 Cover

Bundle MSRP

CA\$20,440

Promotional Code

RL25



- 2 x Rodin RapidCeram A1 300g
- 1 x Rodin Palette 2.0 with Glaze N2-Free

CA\$2,141

Value Added Package

To verify eligibility for the Ackuretta's "Empower Your Workflow with Rodin and SOL Plus" Bundle, including 2 bottles of Rodin RapidCeram A1(CA\$1,372) and 1 x Rodin Palette 2.0 with Glaze N2-Free(CA\$769), customers must register their products through Ackuretta's Concierge Service registration within 30 days of installation and provide proof of purchase. A minimum purchase price of CA\$20,440 is required to be eligible for the promotional goods. Taxes and freight are NOT included. This offer is valid from October 1, 2025 to December 31, 2025. Promotional goods will be sent directly to the customer by Ackuretta within 7-10 days after Concierge Service registration. Ackuretta may modify or cancel this offer at its own discretion without notice. Contact your dealer for more information.

Item Number:

SOL Plus with Concierge Service : RS110010241 | CLEANI with Dual Tanks : RS110010246 | CURIE Plus with N2 Cover : RS110010162



Rodin® SOFT RELINE

Compatible with 3D, PMMA & Acrylic Denture

NEW!

Rodin® Soft Reline features a unique silicone-based formulation that improves the fit and enhances the performance of both removable and fixed full or partial dentures. Tailored for superior bonding and tear strength, it is the ideal reline solution for 3D printed polymer bases, as well as PMMA milled and traditional acrylic materials.

Benefits & Features

- Exceptional Adhesion and Tear Strength
- Ease of Workflow
- Superior Performance
- Esthetic and Patient Comfort

REF 24006 Soft Reliner Kit



Components



VPS Reline Cartridge

Unique silicone-based reline formulation of VPS Reline Cartridge improves the fit and enhances the performance of both removable and fixed full or partial dentures. Its physical properties make it the ideal solution for 3D-printed polymer bases, as well as PMMA milled and traditional acrylic materials.

REF 24007 50mL Reline Silicone Cartridge


REF 240013 Wide Mouth Mixing Tips Refil, 48/pk

VPS Adhesive Primer & VPS Release

VPS Adhesive boasts superior adhesion, thanks to its innovative formulation that penetrates deeply into denture materials to form a resilient bond. Meanwhile, VPS Adhesive Release dissolves tough adhesives without residue, ensuring a comfortable and hassle-free cleaning experience for users.

10mL VPS Adhesive Primer **REF** 24008

10mL VPS Release **REF** 24009




Polishing Brush & Tapered Stone

Specialized polishing tools are included to complement the needs of denture maintenance. The polishing brush assists in refining the reline surface to a flawless finish, while the tapered stone bur shapes materials with precision.

REF 240010 Tapered Cone Finishing Stone

REF 240011 Brown Polishing Wheel

AS SMART AS YOU: UCAN PRINT SPLINT



For clever minds: With UCAN Print Splint by Candulor you get a biomimetic 3D printing resin for the generative production of hard-elastic orthodontic splints featuring a body temperature-induced thermo effect for maximum wearing comfort – and thus a perfect upgrade for your digital workflow.

www.ucan-print.candulor.com

PCU vario 160

Engineered by Dreve



Excellence in Post Curing

Laboratory device on LED basis for curing 3D printed components. Ensuring mechanics and biocompatibility through effective deep curing. The PCU vario offers both the option of a vacuum atmosphere and the option of curing under nitrogen. Biocompatible curing without an inhibition layer can be ensured.

Features

- + Biocompatible curing without inhibition layer
- + Modern, modular system with up to three polymerization chambers
- + Just one control module needed
- + Long-life LEDs
- + Easy to calibrate
- + Very homogeneous intensity
- + Bi-directional illumination from top and bottom
- + 30 storage spaces free to programme
- + Logging and monitoring of process parameters
- + Tray with safe-close function for hermetic sealing of the post curing chamber (vacuum / under nitrogen atmosphere)
- + Unique design: Dreve Red & Stone Grey



Type	Light polymerization unit
Power	1000 W
Weight	45 kg
LED-wave length	405 nm
Light intensity max.	max. 100 mW/cm ² per side
Curing chamber vol.	120 x 480 x 290 mm ³ (H x W x D)
Usable vol.	80 x 420 x 270 mm ³ (H x W x D)
Dimensions of device	240 x 660 x 383 mm ³ (H x W x D)
Environment while curing	Nitrogen, Ambient air, Vacuum





SAVE
\$100

NSK Volvere I7
#Y1002892
SPECIAL PRICE - \$1199.00



\$165
OFF

NSK Ultimate XL Torque
#Y141-466
SPECIAL PRICE - \$1995.00



Neycraft Multizone Duplicating Unit

An innovative and user-friendly unit capable of meeting the changing needs of your lab as workflow increases and decreases.

- Gives user complete temperature control
- Smart capability to work with variable capacities
- Features a safety switch that automatically turns off the mixer if the lid is raised
- Automated power saving on/off feature for an extended motor life
- Comes preprogrammed and ready-to-use

#535631000 3 Gallon \$8995.00



SAVE
\$200

ENBIO S B Class Autoclave
Worlds Fastest Class B Sterilizer

- Simple, fast, small, light and quiet.
- Up to 4 times faster cycle time than the typical autoclave. Easy to use with a small footprint.
- Dry pouches in 15 minutes

#181133460US
SPECIAL PRICE - \$4399.00



SAVE
\$200

Hannox ALGINATE MIXER MX200

Work with perfectly mixed, bubble-free alginate every time with Smile Art Alginate Mixer, a brilliant durable and little benchtop machine.

#ADA004002
SPECIAL PRICE - \$1595.00

Mizzy PIP Pressure Indicating Paste

4+1



Flexos Disposable Saliva Ejectors



4+1

Enamelite[®] Ceramic Spray Glazes



4+1

KeyMill[™] Denture Base Disc



4+1

These offers cannot be combined with any other promotion or discount. Other restrictions may apply. Customer will be invoiced for purchased goods from preferred distributor and no charge goods will ship with the same order. Valid in Canada only.

Promo is valid from 10/1/2025-12/31/2025.



Qu resin

Quick-setting, self-curing denture repair resin based on diacrylate for intraoral and extraoral use.

Refill Cartridge #54001165 - Pink
#54001165 - Dentine



Exakto form

Model resin for accurate reproduction and maximum edge stability.

#52000282- pkg(6) Ivory



bre.Lux Power Unit 2

The bre.lux Power Unit 2 cures all known light curing materials on the market that cure at a wavelength between 370 and 500nm. 14001000

Renfert

100 years making work easy
3+1 year guarantee

Celebrating 100 years of making work easy



SILENT XS

Dust collector for chairside grinding

Item No.: 29221000

\$797.⁵⁰



Dimensions

Height: 6 in.
Length: 7 in.
Width: 9 in.
Weight: 2.5 lbs.



Basic eco

Compact, 2-tank sandblaster

Item No.: 29494025

~~\$1,560.⁰⁰~~

\$1,260.⁰⁰

Renfert



Renfert Firing paste

Firing paste for stabilizing ceramic works

Item No.: 19220304 (4 x 12g)

Renfert's new Firing Paste is the smart choice for dental technicians who demand consistency, clean results, and zero compromise — even with the most delicate ceramics. Usable up to 1,500°C.

\$52.⁹⁵ (4 pk)



Power pillo

Deflasking chisel

Item No.: 50225000

~~\$1,295.⁰⁰~~

\$1,022.⁰⁰



FREE

Flat chisel, wide

With purchase of

Power pillo

Item No.: 50220200

~~\$128.⁹⁵~~





K.2030 Portable BRUSHLESS Rotary Micromotor Kit

The Foredom K.2030 Brushless Portable Micromotor Kit is the upgraded replacement for our K.1030 model. Featuring a brushless motor for longer life, quieter operation, and improved torque, the K.2030 delivers powerful performance in a compact, cordless design.

#K2030

Special Price: \$795.⁰⁰



M.BL Bench Lathe

Powerful, variable speed bench lathe in compact size for use at the workbench.

- Two Year Limited Warranty.
- 1/6 horsepower, high torque motor
- Speed Range: 500-7,000 rpm
- Motor shaft diameter: 5/16"
- Height: 5-1/8", Width: 13" between spindles
- Weight: 7.5 lbs

#MBL

Special Price: \$395.⁰⁰

Aidite

ZIRCONIA PROMOTIONS

10+2



3D PRO ZIR

Buy 10 Get **2 FREE***

*98 X 14mm only

3+1



AIZIR

Buy 3 Get **1 FREE**

*Offer valid through November 30, 2025. Cannot be combined with other promotions. Valid for Aidite Canada customers only.

HUGE

NEW!

4+1



PERFIT Super Light Body

is a high-flow wash material (type 3) that delivers excellent detail reproduction.

It is especially recommended for full denture impressions.

#PER008

Special Price: \$49.⁹⁵

Regular Price: \$7685.⁰⁰

HUGE Multilayer PMMA discs

take advantages of high-quality raw materials and our mature manufacturing technology to realize seamless transition like natural teeth. Multilayer PMMA disc is the most suitable products for aesthetic restorations.



2+1

VITA Removable Prosthetic Deals

Save When You Switch to

VITA



Switch to VITA MFT Teeth from a competitor brand with a \$1500+ NEW stock purchase, get a VITA MFT Living Mould Guide and a **VITA Toothguide® 3D-MASTER** shade guide free*

QTHM2542



Switch to a VITA Premium toothline from a competitor brand with a \$3500+ NEW stock purchase, get a corresponding Living Mould Guide free*

QTHE2542

*one redemption per customer



Buy a **VITA VITACOLL®** Bonding Agent, **OACC2541** get one free

OAKL2543



Buy a **VITA AKZENT® LC** GINGIVA KIT, get a **VITA AKZENT LC GLAZE** free



Buy 4 bottles of **VITA AKZENT LC**, get one free **OAKL2541**

QBAS2541



Buy 10 **VITA VIONIC®** BASE DISCS, get 2 free

QETC2541



Buy 3 **VITA CERAMIC ETCH**, get one free



VITA VIONIC® VIGO Teeth Deals
Buy 20 packs, get a **VITA VIONIC BOND multiDose COMPLETE**
Buy 6 packs, get a **VITA VIONIC BOND Kit**

QVIG2542

QDEN2542

Buy 10 **VITA VIONIC®** DENT DISCS, get a **VITA VIONIC BOND multiDose COMPLETE** free



QFVI2543



Buy a **VITA V60i-Line® PLUS**, get a **VITA AKZENT PLUS MINI KIT** free



QSHD2541



NEW DESIGN!

Buy 2 **VITA Shade Guides**, get one free*

*Free Shade Guide of equal or lesser value



PROMO CODE: Q42025CAN. Valid on purchases made October 1 to December 31, 2025 only. Proof of purchase on one invoice required for freegoods. Freegoods must be equal or lesser value than purchased product. **Invoice must be submitted to promo@vitanorthamerica.com** no later than January 15, 2026. Offers valid in Canada only. Offers may not be combined with any other promotion or discounts. VITA North America reserves the right to discontinue or change the offers at any time without notice. No cash value.

VITA



Scan Here to
Find a Rep
Near You!

EXPERIENCE UNMATCHED PRECISION WITH THE AIDITE EZPRINT-F1 3D PRINTER



**Special
Price**

\$4995.00

Includes:

- Curing & Clean Unit
- 1L Bottle of Model Resin



Ultra high precision

16K ultra-high-resolution LCD printing screen, 14 X 19µm pixel size, restoring design details



Overall uniform lighting

High-power light source system combined with light intensity calibration technology, achieves **95%** light uniformly to control ensuring full-size printing consistency



Efficient production

Large-size platform with **high-speed printing technology** to meet mass production needs




Easy to use

One-touch automatic slicing and printing, **quick-release platform and resin tank**, constant temperature control

EASILY MEET THE NEEDS OF IMPLANTS AND RESTORATIONS

- The crown fits accurately without looseness
- Fit adjustments for analogs in implant models
- 94% overlap rate between dental model and design file*1
- Model deformation rate <3% after 15 days*2



*1: Scanning after the dental model is printed, and comparing with the design model data
*2: Scanning after the dental model is printed, and scanning again after 15 days, compare these two data, the overlap rate value under the tolerance of 0.05mm is ≥97%.